



CLINI*health*[®]

PRODUCT CATALOGUE





DIAGNOSTICS DRUG TESTS & CONSUMABLES

Pregnancy	Page 3
Drug Tests	Page 4 -5
Alcohol Tests	Page 5
Infectious Diseases	Page 6
Urinalysis	Page 7
Lancets and Blades	Page 7 - 8
Intravenous	Page 9
Consumables	Page 10 - 13
Apparatus	Page 13
Ultrasound	Page 14 - 15
Dental Care	Page 15
Bandages	Page 16 - 17
Safety Wear	Page 18
Wound Dressing	Page 19 - 20
Tapes	Page 21
Air Management	Page 22 - 23
Needles and Syringes	Page 24 - 25
Gloves	Page 26 - 27
First Aid Kits	Page 28 - 29
Cotton Products	Page 29
Linen Savers and Diapers	Page 30
Measuring and Checking Apparatus	Page 31 - 36
Packaging and Storage	Page 36

MOBILITY

Bathroom Safety	Page 38 - 41
Rollators	Page 42 - 43
Walkers	Page 44
Wheelchairs	Page 45 - 48
Beds and Mattresses	Page 49 - 54
Stretchers	Page 55 - 57



MOBILITY CONTINUED

Medical Furniture	Page 58 - 60
Trolleys	Page 61 - 64
Crutches + Walking Sticks	Page 65 - 67

ORTHOPAEDICS - ORTHO LIFE

Knee and Leg Support	Page 69 - 71
Back and Shoulder Support	Page 71 - 73
Arm Wrist and Hand Support	Page 73 - 75
Neck and Coccyx Supports	Page 75 - 76
Hernia Belts	Page 76

ORTHOPAEDICS

Compression Stockings	Page 77
Foot and Ankle	Page 77 - 79
Back	Page 80
Hand	Page 80 - 81
Arm	Page 81
Shoulder	Page 82



CLINI *health*
leaders in healthcare



DIAGNOSTICS & CONSUMABLES



Pregnancy



LLOVA005



Easy use
High sensitivity
Early detection
Results in minutes
Assists women in pregnancy planning
ISO and CE approved
FDA Approved
Comes with monthly planner
5 tests per box

Barcode:700371643867
Nappi Code:184627001

**Ovulation Test
5 pack**

CHPTD50



50 Individually Packed Devices.
Colloidal Gold Technology.
High Accuracy 99%, Packed in 50s.
One Step Easy Use.
CE TUV and FDA Certified Approved.
Pregnancy tests are affordable and highly accurate means of testing for early pregnancy.

50 x Individually wrapped tests

Barcode:6009701458268
Nappi Code:1004137001

Pregnancy Device - individually packed retail - 50s

EPT001



Colloidal Gold Technology.
High Accuracy 99%.
One Step Easy Use.
CE TUV and FDA Certified Approved.
Pregnancy tests are affordable and highly accurate means of testing for early pregnancy.

1 individually wrapped test

Barcode:700371643843
Nappi Code:184625001

Pregnancy Midstream Test - Singles

EPTM25



25 Individually Packed Tests.
Colloidal Gold Technology.
High Accuracy 99%, Packed in 25s.
One Step Easy Use.
CE TUV and FDA Certified Approved.
Pregnancy tests are affordable and highly accurate means of testing for early pregnancy.

Barcode:700371643850
Nappi Code:184625002

Pregnancy Midstream test Shipper - 25s

BOVPREGTEST1




Test time: 20 min, room temperature
Sample type: Serum, heparin sodium plasma, or EDTA plasma of cows
Early diagnosis: Cows can be tested after 28 days post-breeding.

Barcode: 0710579687365
Nappi Code:

Cow Pregnancy Tests Bovine – Singles

EPT002



Double Pack of tests with free cup
Easy use
High Sensitivity
Early Detection
Results in Minutes
For Self Test
Manufactured under ISO 13485
Each Box contains 2 Pregnancy Tests

Barcode: 700371643836
Nappi Code: 184626001

Pregnancy Test Double Pack

CHPTS25




Colloidal Gold Technology, High Accuracy 99%, Packed in 25s,
One Step Easy Use, CE TUV and FDA Certified Approved,
Pregnancy tests are affordable and highly accurate means of testing for early pregnancy.

Barcode:6009701458275
Nappi Code:1004139001

Pregnancy Test Strip 25s

CHPTS100




Colloidal Gold Technology, High Accuracy 99%, Packed in 100s,
One Step Easy Use, CE TUV and FDA Certified Approved,
Pregnancy tests are affordable and highly accurate means of testing for early pregnancy.

Barcode:6009701458282
Nappi Code:1004139002

Pregnancy Test Strip 100s

Drug Tests



D5CH1


High Sensitivity
Advanced Quality

Tests for:
Cocaine / COC
Amphetamine / AMP
Dagga-Cannabis / THC
Opiates – Morphine – Heroin / MOP
Methamphetamine / Tik / Crystal Meth / MET

Available in packs of 10s

Barcode:6009879738292
Nappi Code:210491001

Drug Test 5 Single



D6CH1


High Sensitivity
Advanced Quality

Tests For: Cocaine / COC
Amphetamine / AMP
Dagga-Cannabis / THC
Opiates – Morphine – Heroin / MOP
Methamphetamine / Tik / Crystal Meth / MET
Benzodiazepine – Valium and other sedatives / BZO

Available in packs of 10s

Barcode:6009879738308
Nappi Code:210566001

Drug Tests 6 Single



D17CH1

Multiple Drug Tests in One Go

Fast results

Easy visually interpretation

Simple operation, no equipment required

High accuracy

Barcode:6009701456226
Nappi Code:1143550001

Drug Test 17 Clinihealth



D12CH1

High Sensitivity
Advanced Quality
ISO Approved
SAMHSA Standard
Results in minutes
Convenient Saliva Type

Barcode:6009701456332
Nappi Code:N/A

Drug Test 12 Saliva



DRUGCUP6

Clinihealth drug test 6 panel cup.
Tests For: Cocaine / COC
Amphetamine / AMP
Dagga-Cannabis / THC
Opiates – Morphine – Heroin / MOP
Methamphetamine / Tik / Crystal Meth / MET
Benzodiazepine – Valium and other sedatives / BZO

Barcode:6009701452280
Nappi Code:1087723001

Drug test 6 Panel Cup



DRUGCUP12

High Sensitivity, Advanced Quality
Results in minutes, Tests for 12 different drug types

AMP – Amphetamine	MTD – Methadone
BAR – Barbiturates	MOP – Morphine
BZO – Benzodiazepines	OPI – Opiates
COC – Cocaine	THC – Cannabis
MDMA – MDMA/Ecstasy	TRA – Tramadol
MET – Methamphetamine	
TML – Tamsulosin Hydrochloride	

Barcode:6009701451009
Nappi Code:1087715001

Drug test 12 Panel Cup



DRUGCUP17

CliniHealth Drug Test 17 panel CUP tests for 17 different types of abused drugs
High Sensitivity, Advanced Quality.
ISO Approved, CE Approved, FDA Approved, SAMHSA Standards.

Results in minutes.

Barcode:6009701456233
Nappi Code:N/A

Drug test 17 Panel Cup

Drug Tests

CHTHC50



High Sensitivity, Advanced Quality, ISO Approved, CE Approved, FDA Approved, SAMHSA Standards, Results in minutes, Tests for THC (Dagga)

Barcode :6009701450705
Nappi Code:1087727001

THC Test Clinihealth - 50s

DRUG6SALCH



High Sensitivity
Advanced Quality
ISO Approved
CE Approved
SAMHSA Standard
Results in minutes
Tests for 6 different types of commonly abused types
Convenient Saliva Type

Barcode :6009701452648
Nappi Code:1129554001

Saliva Drug Test 6 Panel

D8CH1



High Sensitivity, Advanced Quality, ISO Approved, CE Approved, FDA Approved, SAMHSA Standards, Results in minutes,

Available in packs of 20s

Opiates - OPI
Methamphetamine / Tik / Crystal Meth -MET
Dagga-Cannabis - THC
Methadone - MTD
Ecstasy -MDMA
Benzodiazepine – Valium and other sedatives -BZO
Cocaine - COC
Opiates – Morphine – Heroin - MOP

Barcode :6009701450965
Nappi Code:1129549002

Drug Test 8 Panel

DRUGCUP12A



High Sensitivity, Advanced Quality
Results in minutes, Tests for 12 different drug types

AMP – Amphetamine	MTD – Methadone
BAR – Barbiturates	MOP – Morphine
BZO – Benzodiazepines	OPI – Opiates
COC – Cocaine	THC – Cannabis
MDMA – MDMA/Ecstasy	ALC - Alcohol
MET – Methamphetamine	
TML – Tamsulosin Hydrochloride	

Barcode: N/A
Nappi Code:

Drug test 12 Panel Cup with Alcohol

Alcohol Tests



ALURINCHI

The Alcohol Rapid Test is a rapid, highly sensitive method to detect the presence of alcohol in specimen and provide an approximation of relative alcohol concentration.

Fast results

Easy visually interpretation

Simple operation, no equipment required. **Available in packs of 40s**

High accuracy

Barcode:6009701456240

Nappi Code:1143424001

Alcohol Rapid Test Urine Device Singles



HALCBRT

The Breath Alcohol Test is a preventive device intended to give you an indication of your blood alcohol level.

Tests 4 different levels of alcohol in breath : 0,00 %, 0,02% 0,05% 0,08%.

Barcode: 6009701459371

Nappi Code:

Alcohol Breath Test - Singles

The Alcohol Rapid Test is a rapid, highly sensitive method to detect the presence of alcohol in specimen and provide an approximation of relative alcohol concentration.



OST25CH

Barcode: 0710579687372

Nappi Code:

One Step Ethyl ETG Alcohol Test Urine - 25s



ALORALCHI

The Alcohol Rapid Test is a rapid, highly sensitive method to detect the presence of alcohol in specimen and provide an approximation of relative alcohol concentration.

Fast results

Easy visually interpretation

Simple operation, no equipment required. **Available in packs of 25s**

Simple operation, no equipment required

Nappi Code:1143426001

Alcohol Rapid Test Oral Device Singles

Infectious Disease Tests



CHHIV1

HIV tests are the most commonly used means of controlling and managing the spread of HIV. The HIV test most commonly used is the rapid HIV TEST kit. The best way to test for HIV is a rapid test in order to combat and control the HIV epidemic. These HIV tests look for the antibody that the immune system produces.

Barcode :6009879738377
Nappi Code:274350001

ing test
CPV 380 Anti-HIV
CPX 070 Anti-HIV

Confirmatory test
CPX 025 Anti-HIV (V

HIV-Test

Hiv Kit 1&2 Comp Clinihealth (incl swab /lancet/plaster)



MAND25CH


The Clinihealth Mandrax MQL tests for Methaqualone (Mandrax). High Sensitivity. Advanced Quality. ISO Approved, CE Approved, FDA Approved, SAMHSA Standards.

Results in minutes

Barcode :6009701450033
Nappi Code:1116935001



Mandrax Test 25s (MQL)




CHMAL1

The Clinihealth Malaria Test Device complete single, contains swab, lancet and plaster. The rapid test simultaneously tests for the fatal strains of PF as well as PV – Plasmodium faciparum and Plasmodium vivax.

Results in minutes.
Trusted in 64 Countries.
High Sensitivity.
Advanced Quality.
High Specificity.
Available in packs of 25s

Barcode :6009879738285
Nappi Code:196710001



Malaria Test Device Single



CHSYPC25

Clinihealth Syphilis Test. Clear preparation procedures. Simple and safe specimen collection directions. Comprehensive package of required materials and components. Rapid and accurate test results within 10-20 minutes. High Sensitivity. Advanced Quality. Reliable and inexpensive. Convenient and easy-to-understand directions.

Barcode :6009701450736
Nappi Code:1087766001



Syphilis Test 25s Available in singles and 25s

Urinalysis



URICH6P

High Sensitivity.
Advanced Quality.
ISO Approved.
CE Approved.
FDA Approved.
SAMHSA Standards.
6 Parameter Professional.
GP Urinalysis Multisticks.
Urine Strip Test Stick.
Strips for Leucocytes, Nitrites, Urobilogen, Protein, Glucose.
Pack of 100 strips

Barcode :6009701458299
Nappi Code:1004143001

Urinalysis Combi 6




URI10CH

High Sensitivity.
Advanced Quality.
ISO Approved.
CE Approved.
FDA Approved.
SAMHSA Standards.
10 Parameter Professional.
GP Urinalysis Multisticks.
Urine Strip Test Stick Strips for Blood, Billirubin, Urobilogen, Ketone, Protein, Nitrite, Glucose, PH, Specific Gravity, Leucocytes.
Pack of 100 strips

Barcode : 6009879738407
Nappi Code: 274359001

Urinalysis Combi 10



URINALYSERW200

Test Items: Leucocytes, Nitrite, urobilinogen, protein, pH, blood, specific gravity, VC, ketone, bilirubin, glucose, micro albumin- Strips: URS-10T, H11-MA, H12-MA- 60 strips/h for single test, 120 strips/h for consecutive test- Memory: 1000 results, including time and patient ID- Built in printer

Barcode:6009701450019
Nappi Code:N/A

Urine Analyser W200B CliniHealth



Lancets and Blades



CHLAN10023

Easy Comfort
Gamma sterilized
Advanced Needle Tip

- Advanced Design
- Less Needle Failure
- Ergonomic Design
- Gamma Sterilized

Barcode :8906061780035
Nappi Code:274352-001

Safety Lancets CliniHealth (23g)



CHLAN10028

Easy Comfort
Gamma sterilized
Advanced Needle Tip

- Advanced Design
- Less Needle Failure
- Ergonomic Design
- Gamma Sterilized

Barcode :6009701450798
Nappi Code:1087731001

Safety Lancets CliniHealth (28g)



LANBLD200

Blood Lancets
200s

- Sterile, disposable, single use lancets
- Designed for easy puncturing to minimize trauma
- Tapered stainless steel point
- Consistent depth penetration
- Ideal for finger sticking or finger pricking
- Designed for easy puncturing to minimize trauma STEEL

Barcode :6009701453973
Nappi Code:1135261001

Stainless Steel Blood Lancets -200s



CHLAN10021

Easy Comfort
Gamma sterilized
Advanced Needle Tip

- Advanced Design
- Less Needle Failure
- Ergonomic Design
- Gamma Sterilized

Barcode :6009701451450
Nappi Code:1108754001

Safety Lancets CliniHealth (21g)

Lancets and Blades



BLADSURG10



BLADSURG12



BLADSURG11



BLADSURG15

Disposable Carbon Steel scalpel blades
High quality Carbon steel blade, Convex cutting surface.
Non-attached, replaceable blades, easy to insert into handle.
These blades are individually packed in peel-open pouches.

No 10
Barcode:6009701454604
Nappi Code:1135272001

No 11
Barcode:6009701454611
Nappi Code:1135273001

No 12
Barcode:6009701454642
Nappi Code:1135274001

No 15
Barcode:6009701454666
Nappi Code:1135275001

Carbon Steel Surgical Blades - 10/11/12/15



TLAN30G

TWIST TOP: The protective twist top allows for ease of use and a secure way of disposal.

After testing, reattach the cap and dispose of properly to prevent dangers posed by the sharp needle and the possible spread of disease-causing agents.

These lancets are for single use only.

Barcode :6009701451054
Nappi Code:1135262001

Twist Top Lancet



Made of high grade non rusted stainless steel
Can be autoclave
Length- No.3-120mm
Length- No.4-130mm
Quality Products.

SCHDLCH3
No 3

Barcode :6009701457445
Nappi Code:1185610001

SCHDLCH4
No 4

Barcode :6009701457452
Nappi Code:1185611001

Scalpel Handle Stainless Steel No 3 & 4



Scalpel blades can has been used: in seedling inoculation to dissect adrenal gland for adrenalectomy

To cut the polydimethylsiloxane (PDMS) layers for fabrication of the microfluidic device

No 10 - SURBLDPLH10
Barcode:6009701457384
Nappi Code:1155590001

No 20 - SURBLDPLH20
Barcode:6009701457414
Nappi Code:1155577001

No 11 - SURBLDPLH11
Barcode:6009701457391
Nappi Code:1155592001

No 21 - SURBLDPLH21
Barcode:6009701457421
Nappi Code:1155581001

No 12 - SURBLDPLH12
Barcode:6009701457001
Nappi Code:1155574001

No 23 - SURBLDPLH23
Barcode:6009701457438
Nappi Code:1155596001

No 15 - SURBLDPLH15
Barcode:6009701457407
Nappi Code:1155574001

Surgical Blade Stainless Steel with Plastic Handle - 10/11/12/15/20/21/23

Intravenous



Single Use
Sterilized by EO
Pyrogen Free
Disposable Sterile

100 x Sterile Radiopaque Catheters

- Store in a cool dry place and away from direct sunlight
- Do not use if packaging is damaged
- Discard after single use
- No Port and Wings

18G - IVC18G1345CH

Barcode:6009701457155

Nappi Code:1143439001

20G - IVC20G1132CH

Barcode:6009701457131

Nappi Code:1143437001

22G - IVC22G0925CH

Barcode:6009701457117

Nappi Code:1143436001

24G - IVC24G0719CH

Barcode:6009701455687

Nappi Code:1143434001

IV Catheter without port and wings - 100s



IVSETAD201X50

20 Drop Adult
2 piece rotary luer lock
Contains 1 latex free Y-Site
Sterilized by E.O. Gas
Non Pyrogenic
Latex Free
Individually wrapped in a peel-open blister pouch
For gravity use only

Barcode :6009701455946

Nappi Code:1143447001

IV Admin Set Adult 20 Drops 1Y Site - 50s



IVSETAD601X50

60 Drop Peads
2 piece rotary luer lock
Contains 1 latex free Y-Site
Sterilized by E.O. Gas
Non Pyrogenic
Latex Free
Individually wrapped in a peel-open blister pouch
For gravity use only

Barcode :6009701450026

Nappi Code: N/A

IV Admin Set Pediatric 60 Drops 1Y Site

Single Use
Sterilized by EO
Pyrogen Free
Disposable Sterile

100 x Sterile Radiopaque Catheters

- Store in a cool dry place and away from direct sunlight
- Do not use if packaging is damaged
- Discard after single use
- With Port and Wings



18G - IVCPW18G1345CH

Barcode:6009701457179

Nappi Code:1143446001

20G - IVCPW20G1132CH

Barcode:6009701457148

Nappi Code:1143444001

22G - IVCPW22G0925CH

Barcode:6009701457124

Nappi Code:1143442001

24G - IVCPW24G0719CH

Barcode:6009701455700

Nappi Code:1143440001

IV Catheter with port and wings - 100s

Intravenous

FOLCT2W10CH



Foley catheter used for short or long term urine drainage.

Made up of latex bonded with silicone elastomer which eliminates the risk of encrustation.

Smooth surface facilitates non-traumatic incubation.

Barcode :0710579688355
Nappi Code:

Foley Catheter Latex Silicone Two Way 10FR - Singles

FOLCT2W12CH



Foley catheter used for short or long term urine drainage.

Made up of latex bonded with silicone elastomer which eliminates the risk of encrustation.

Smooth surface facilitates non-traumatic incubation.

Barcode :0710579688362
Nappi Code:N/A

Foley Catheter Latex Silicone Two Way 12FR - Singles

FOLCT2W16CH



Foley catheter used for short or long term urine drainage.

Made up of latex bonded with silicone elastomer which eliminates the risk of encrustation.

Smooth surface facilitates non-traumatic incubation.

Barcode :0710579688386
Nappi Code:

Foley Catheter Latex Silicone Two Way 16FR - Singles

FOLCT2W14CH



Foley catheter used for short or long term urine drainage.

Made up of latex bonded with silicone elastomer which eliminates the risk of encrustation.

Smooth surface facilitates non-traumatic incubation.

Barcode :0710579688379
Nappi Code:

Foley Catheter Latex Silicone Two Way 14FR - Singles

FOLCT2W20CH



Foley catheter used for short or long term urine drainage.

Made up of latex bonded with silicone elastomer which eliminates the risk of encrustation.

Smooth surface facilitates non-traumatic incubation.

Barcode : 0710579688409
Nappi Code:

Foley Catheter Latex Silicone Two Way 20FR - Singles

FOLCT2W18CH



Foley catheter used for short or long term urine drainage.

Made up of latex bonded with silicone elastomer which eliminates the risk of encrustation.

Smooth surface facilitates non-traumatic incubation.

Barcode :0710579688393
Nappi Code:

Foley Catheter Latex Silicone Two Way 18FR - Singles

FOLCT2W22CH



Foley catheter used for short or long term urine drainage.

Made up of latex bonded with silicone elastomer which eliminates the risk of encrustation.

Smooth surface facilitates non-traumatic incubation.

Barcode :0710579688416
Nappi Code:

Foley Catheter Latex Silicone Two Way 22FR - Singles

Consumables

CHTD100

- Wooden tongue depressors.
- Packed in 100's,
- Ergonomic design,
- Supreme A grade wood
- Easy Grip.



Barcode:6009701452136
Nappi Code:N/A

Tongue Depressor

CHAW200

Alcohol swabs are cotton pads coated in pure isopropyl alcohol, sterile and individually packed. They are generally used to sterilize a site prior to injection and application. Ideal in preventing bacteria and germs from entering and infecting wound sites resulting in infection. Each pouch is sealed in a sterile aluminium pouch to prevent drying out- other uses include cleaning small cuts scrapes and blisters.



Barcode: 6009701452228
Nappi Code: 274348-001

Alcohol Swabs - 200s

Single Use
Sterile
Individually Packed
Latex Free



VAGSPC100CHS
Barcode:N/A
Nappi Code:710579687242

VAGSPC100CHM
Barcode:6009701455809
Nappi Code:1147248001

VAGSPC100CHL
Barcode: 6009701458183
Nappi Code: 1166346001

Vaginal Speculum with Screw Type

URNBG2000

Urine Drainage Bag with Screw Valve (twist valve) is made of polyvinyl chloride (PVC) with medical grade. T-Tap cross outlet valve Non-return valve 2L Comes in a Sterile pouch pack - Pack of 10s



Barcode:6009701454024
Nappi Code:1135260001

Urine Bag- 2000ml - 10s

KIDDISHSS1

25x10x5cm
Stainless Steel



Barcode:6009701457469
Nappi Code:N/A

Kidney Dish Stainless Steel - 25cmx10cmx5cm

TOUQTCH1

Tourniquet Large With Buckle




Barcode:6009701452617
Nappi Code:

Tourniquet with Buckle

CHENBI

Endometrial investigation:
- Detection of luteal phase insufficiency,
- Monitoring of hormonal treatment,
- Abnormal bleeding,
- Bacteriological cultures to identify pathogens.




Barcode: 6009701459401
Nappi Code:N/A

Endometrial Biopsy Curetter - 20s

Important way to reduce patient safety risks
Comfortable to wear
Premium quality

Soft vinyl with adjustable snaps
Hospital use
Includes patient information



Pink- Adult Product Code: IDBND100P Barcode:6009701453768 Nappi Code: 1135259001	White - Neonatal Product Code: IDBND100CW Barcode:6009701450989 Nappi Code: 1166353001	Pink - Neonatal Product Code: IDBND100NP Barcode:6009701458732 Nappi Code: 1166317001	Pink- Child Product Code: IDBND100CP Barcode:6009701459756 Nappi Code: 1166349001
Blue - Neonatal Product Code: IDBND100NB Barcode:6009701458749 Nappi Code: 1166317001	Blue - Adult Product Code: IDBND100B Barcode:6009701453744 Nappi Code: 1135258001	Yellow - Child Product Code: IDBND100CY Barcode:6009701458763 Nappi Code: 1166351001	White - Adult Pink: IDBND100W Barcode:6009701458725 Nappi Code: 1166347001

Identification Bands - 100s

Consumables

SKSTCH01

Skin Stapler 6.5mm wide 35 Staples in a Box | Disposable



Barcode: 6009701458893
Nappi Code: 1185612001

Skin Stapler

SKSTRCH02

Highly Efficient. Avoid Injury To Epidermal Tissue. Simple and Reasonable Design. Single Use Disposable.



Barcode: 6009701458909
Nappi Code: 1185613001

Skin Stapler Remover

BIOHAZBCH

They are used to dispose of used micro tips, tubes and other plastic products and contaminants. The bags are designed to withstand bio hazardous waste and are clearly marked with a 'Bio Hazard' print. The bags are made from autoclavable polypropylene. 35cm x100cm - 20s

Available also in Autoclavable format **BIOHAZACH**



Barcode: 6009701450002
Nappi Code: N/A

Bio Hazard Bag - Plastic - 35cm x100cm - 20s

5LSHCCH

Made of Thick and durable plastic Leak proof and puncture proof Warning signs that indicate that container has biohazard and sharp objects Top handle for easy transportation Rotating top cap to help secure needles.

Specifications:
Dimension: 26 x 21 x 18 cm



Barcode: 6009701458930
Nappi Code: N/A

5L Sharp Container

MEDRSCH

Stainless steel blade Sturdy plastic handle Blunt tip with snip protection Convenient angled design Excellent for cutting dressing and bandages



Barcode: 6009701458916

Medical Rescue Scissors 19cm

ICBBRCH

2cm (D) 9.5cm (W) 14cm (H) 200gm ice pack Simply place in the freezer before use Leak proof cap to prevent spillage



Barcode: 6009701458923

Ice Brick Blue 200g

CHVMT

Plastic material Product size: 21x15x9cm



Barcode: 6009701458916

Vomit Bowl

URB1200

Blue 1200ml



Barcode: 6009701459388
Nappi Code: 1224038001

Urinal Blue 1200ml

Consumables

DISSFOT14

Straight, serrated forceps is excellent for removing surgical dressings or for other detail work requiring a precision tweezer.



Barcode: 6009701459661

Dissection Forceps with Teeth 14cm

DISSFOW14

Straight, forceps is excellent for removing surgical dressings or for other detail work requiring a precision tweezer.



Barcode: 6009701459678

Dissection Forceps without Teeth 14cm

HEMFST14

Kelly Straight 14cm stainless steel Hemostatic Forceps




Barcode: 6009701459685

Hemostatic Forceps Kelly Straight 14cm

SCMSR14

Mayo scissors have semi-blunt ends, a feature that distinguishes them from most other surgical scissors. Metzenbaum scissors are, however, similar in that the ends are the same on both scissor halves, but its handles are longer and its middle section is slightly narrower.



Barcode: 6009701459692

Scissors Mayo Straight 14cm

SCMCR14

Mayo scissors have semi-blunt ends, a feature that distinguishes them from most other surgical scissors. Metzenbaum scissors are, however, similar in that the ends are the same on both scissor halves, but its handles are longer and its middle section is slightly narrower.




Barcode: 6009701459708

Scissors Mayo Curve 14cm

HEMFST18

Kelly Hemostatic Forceps are primarily used for clamping large blood vessels or manipulating heavy tissue. Kelly Forceps may also be used for soft tissue dissection




Barcode: 6009701459715

Hemostatic Forceps Kelly Straight 18cm

20LSHCCH

Medical waste bin, 10 litre capacity for disposing of bio-hazardous needles & sharps.




Barcode: 6009701450996

10L Sharp Container

Consumables

URNEBAGPTB



Skin Stapler 6.5mm wide 35 Staples in a Box | Disposable

Barcode: 0710579686962
Nappi Code: N/A

Portable Urine System

CHKB1



Stainless Steel 201 Frame 304ss Bowl
Bowl Diameter: 28cm

Barcode: 0710579686979
Nappi Code: N/A

Kick Bowl

CHKB2



Stainless Steel 201 Frame 304ss Bucket
Bucket Diameter: 30cm

Barcode: 0710579686986
Nappi Code: N/A

Kick Bucket

CHBPCV



304SS With Cover

Barcode: 0710579686993
Nappi Code: N/A

Bedpan with Cover

CHBPPL



Bedpan

Barcode: 0710579687389
Nappi Code: N/A

Plastic Bedpan

CHKBT



Forceps Plastic are designed for a comfortable and secure grip. Their ergonomic design and plastic construction make them durable and lightweight. Perfect for work in many medical and laboratory applications.

Barcode: 6009701459562
Nappi Code: N/A

Plastic Forceps - 100s

CHBPPL



Our range of disposable bio-hazard bags are bright red in colour and easy to open. They are used to dispose of used micro tips, tubes and other plastic products and contaminants. The bags are designed to withstand bio hazardous waste and are clearly marked with a 'Bio Hazard' print.

Size: 210x310x0.3mm
Colour: bright red

Barcode: 6009701459982
Nappi Code: N/A

Bio Hazard Bag Autoclavable - Plastic - 210x310x0.3mm - 20s

Consumables


CHDS201 Dressing Bowl



Barcode: 0710579688423
Nappi Code: N/A

Dressing Bowl 201SS

CHDS304 Dressing Bowl



Barcode: 0710579688430
Nappi Code: N/A

Dressing Bowl 304SS

CHBT201 Bend Tray 201SS



Barcode: 0710579688447
Nappi Code: N/A

Bend Tray 201SS


CHSST201 Stainless Steel Tray 201



Barcode: 0710579688461
Nappi Code: N/A

Stainless Steel Tray 201

CHBT304 Bend Tray 304SS



Barcode: 0710579688454
Nappi Code: N/A

Bend Tray 304SS

CHSST304 Stainless Steel Tray 304



Barcode: 0710579688478
Nappi Code: N/A

Stainless Steel Tray 304

CHKB1 Stainless Steel Box Disinfection Case Storage Tank for Lab Veterinary Dental Business & Industrial



Available in 17cm / 20cm / 23cm / 26cm / 30cm

Barcode: Available on Request
Nappi Code: N/A

Storage Tank

CHTT304 Stainless Steel Antiodophor Medical Instrument Treatment Tray



Barcode: 0710579688539
Nappi Code: N/A

Treatment Tray 304

Consumables



CHTB304

Tweezers Bucket

Barcode: 0710579688546
Nappi Code: N/A

Tweezers Bucket



Convenient to protect casts or bandages from the water while showering or taking part in water activities
Long Arm – 630mmx 430mm – 180mm ring diameter
NOTICE – Plastic slippery when wet, be careful of wet surfaces or protectors.
Plastic can cause suffocation, keep out of children's reach unless under adults supervision.
Do not expose the protector to heat or flames

<p>Size/Product: Cast Cover Arm - Adult Product Code: CHCVARA Barcode: 0710579687921 Nappi Code: N/A</p>	<p>Size/Product: Cast Cover Arm - Child Product Code: CHCVARC Barcode: 0710579687938 Nappi Code: N/A</p>
--	--

Cast Cover

Apparatus

ELECPNPCCH



Electro Coagulation Pen
Cutta Cutter Handle x 1 Pc
Charger Base
Pen Tips x 4

TIPS CAN BE BOUGHT SEPARATELY

Barcode:6009701458589
Nappi Code:N/A

Cautery Pen - Electro Coagulation Pen

ENTDIAWCH



Technical Configuration of Wall-mounted ENT Diagnostic Set

- Wall-mounted Type
- Ophthalmoscope
- Otoscope
- Nasal speculum
- Electronic thermometer (Forehead and ear type)
- Blood Pressure Monitor
- Storage box
- Hanging plate
- Optional parts: tongue depressor, finger clip type pulse oximeter, blood-glucose meter

Barcode:6009701458954
Nappi Code:N/A

Wall Mounted ENT Diagnostic Set

ENTDIAGCH



Chrome-plated bent arm illuminator for lighting the oral cavity and pharyngeal area
Two chrome-plated laryngeal mirrors
Ophthalmoscope with lens wheel and correcting lens diopters ranging from -20 to + 20
Chrome-plated nasal speculum
Tongue Blade Holder
Handle is operated by two c batteries <not included>
Quality plastic box packing for extra protection.

Barcode:6009701458954
Nappi Code:N/A

ENT Diagnostic Set

CHDIOPTA



Direct ophthalmoscope a classic medical device for ophthalmology adopts a professional lens system.

Barcode: 0710579687006
Nappi Code:N/A

Diagnostic Ophthalmoscope

ENTDIAGCH



1 Otoscope
1 Ophthalmoscope

Barcode: 0710579687013
Nappi Code:N/A

Otoscope Black

Ultrasound

ULTGBG5LCH



- Viscous Gel for general
- Ultrasonography applications
- Use for ECG and ultrasound
- Water Soluble
- Clear

Barcode:6009701458947
Nappi Code: N/A

Ultrasound Gel Bag Clear 5L

Ultrasound

ULTGEL500CH



- Viscous Gel for general use
- Ultrasonography applications
- Use for ECG and ultrasound
- Water Soluble
- Clear
- Easy to use Bottle


Barcode:6009701456912
Nappi Code: N/A

UltraSound Gel 500ml



Dental Care





TOOTBRNY12CH

Various Colours 12s

Tear off sections with individual barcodes

Hanging Card

Barcode:6009701455885
Nappi Code:N/A

Toothbrush - Various Colours 12s

Our premium denture box helps maintain a high level of denture hygiene by providing a clean storage space for your dentures or fake teeth, retainers, or sports or athletic mouth guards. Our strong denture case also makes it convenient to carry dentures anywhere and prevents accidental handling damages.

Black - DENBLKCH
Barcode:6009701457490
Nappi Code:N/A

Blue - DENBLUCH
Barcode:6009701457476
Nappi Code:N/A

Purple - DENPURCH
Barcode:6009701457513
Nappi Code:N/A

Green - DENGGRNCH
Barcode:6009701457483
Nappi Code:N/A



Denture Box

Bandages



This type of bandage is a stretchable bandage and elasticised. Crepe bandages is used to treat muscle sprains and injuries on bones and tendons. By correct application it assists in reducing pain and inflammation by means of evenly distributed pressure. This creates faster recovery and ease of comfort during recovery. Due to its heavier weight it can be used on ankles, knees and many other areas requiring heavy use.

- Packed individually with elastic fabric clip

Crepe Bandage 50mm
x 4.5m - CREBND50
Barcode :6009701452310
Nappi Code:1006027001

Crepe Bandage 100mm
x 4.5m - CREBND100
Barcode :6009701452341
Nappi Code:1006020001

Crepe Bandage 150mm
x 4.5m - CREBND150
Barcode :6009701452334
Nappi Code:1006025001

Crepe Bandage 75mm
x 4.5m - CREBND75
Barcode :6009701452327
Nappi Code:1006026001

Crepe Bandages



A woven, stretchy bandage, for medical use
Each roll individually wrapped in a coloured pouch
For mild compression

CNFBND5X45
Conforming Bandage - 5cmx4.5m
Barcode :6009701454512
Nappi Code:1135263001

CNFBND7545
Conforming Bandage - 7.5cmx4.5m
Barcode : 6009701454536
Nappi Code:1135264001

CNFBND1045
Conforming Bandage - 10cmx4.5m
Barcode : 6009701454550
Nappi Code:1135265001

CNFBND1545
Conforming Bandage - 15cmx4.5m
Barcode :6009701454574
Nappi Code:1135266001

Conforming Bandages



Orthopaedic Padding Bandage, For Medical
Orthopaedic padding is the wide-used supporting material for surgery and swaddling of plaster of paris bandage.
Made of natural or synthetic non-woven fabrics.

ORPAD75X3CH
Orthopaedic Padding 7.5x3m - 12s
Barcode :710579687198
Nappi Code:N/A

ORPAD10X3CH
Orthopaedic Padding 10x3m - 12s
Barcode : 710579687228
Nappi Code:N/A

ORPAD15X3CH
Orthopaedic Padding 15x3m - 12s
Barcode : 710579687235
Nappi Code:N/A

Orthopaedic Bandages



Stretchable extra firm fabric supports and moulds around joints that require support
Ideal for joints that are aching and need to be pampered or reinforced during exercise
Quick and easy to apply – sticks only to itself
Does not stick to hair or skin
Non-slip support – no pins or clips needed
Lightweight and breathable
Multi-use
Ultra flexible providing maximum comfort, maintaining its shape and will not loosen or shift
Colour: natural

Self Adherent Cohesive
2.5cm x 4.5m - SABND25
Barcode :6009701453867
Nappi Code:1135267001

Self Adherent Cohesive
5cm x 4.5m - SABND50
Barcode :6009701453874
Nappi Code:1135268001

Self Adherent Cohesive
7.5cm x 4.5m - SABND75
Barcode :6009701453881
Nappi Code:1135269001

Self Adherent Cohesive
10cm x 4.5m - SABND100
Barcode :6009701453898
Nappi Code:1135270001

Bandage Self Adherent Cohesive

Bandages



Plaster cloth is one of the most often used sculpting products and is essentially the same product used by doctors to make plaster casts. The product is made of pre-shrunk cotton gauze fabric infused with Plaster of Paris. Strong and fast setting. Setting time of 3-4 minutes. Creates strong cast of immobilizing broken and fractured limbs. Excellent conformability and low plaster loss. Accurate and easy moulding with high quality creamy textured Plaster of Paris. Each roll individually sealed in a special moisture resistant flow wrapper for easy storage.

Plaster of Paris Bandage
5 x 3M - POPBND50
Barcode :6009701453904
Nappi Code:1135336001

Plaster of Paris Bandage
7,5 x 3M - POPBND75
Barcode :6009701453911
Nappi Code:1135337001

Plaster of Paris Bandage 10 x 3M
POPBND100
Barcode :6009701453928
Nappi Code:1135250001

Plaster of Paris Bandage
15x 3M - POPBND150
Barcode :6009701453935
Nappi Code:1135251001

Plaster of Paris Bandage 20x 3M
POPBND200
Barcode :6009701453942
Nappi Code:1135252001

Plaster of Paris Bandages



NWTRBNDGE1CH

Multi Purpose - Triangular bandages are commonly used to compress various injuries and used as arm sling, cravat bandage, cover for head dressing, support sprains and/or broken bones. 100% Cotton - Made from top quality 100% non-woven cotton for maximum comfort. Individually Packed - Each arm sling comes in a individually zip closed poly bag making sure the bandages get to you clean and safe to use. Latex Free

Barcode : 6009701455816
Nappi Code:1143551001

Bandage Non Woven Triangular



FADN37575
FADN41012
FADN51216
FADN61517

First Aid Dressing is a hygienic pressure dressing for wounds to minimize bleeding and help keep cuts, scrapes and other injuries clean and dry.

First Aid Dressing No 3
Pad Size 75mmx75mm
FADN37575
Barcode :6009701458312
Nappi Code:1166249001

First Aid Dressing No 5
Pad Size 120mmx160mm
FADN51216
Barcode :6009701458336
Nappi Code:1166251001

First Aid Dressing No 4
Pad Size 100mmx120mm
FADN41012
Barcode :6009701458329
Nappi Code:1166250001

First Aid Dressing No 6
Pad Size 150mmx170mm
FADN61517
Barcode :6009701458343
Nappi Code:1166252001

First Aid Dressing



Flexible and stretchy to support injuries but allow movement. Strong adhesive can be relied upon in training and competition. Soft, feathered edges for comfortable strapping. Easy to unwind for quick and simple application. Centre guideline for even overlap. Each roll individually wrapped for convenience and hygiene. A yellow/red center line on the back of the cloth acts as an overlap guide for uniform application and pressure. The dry feathered edge provides for a neat, comfortable application.

Elastic Adhesive Bandage
2.5 x 3M - EABBND25X3
Barcode :6009701456424
Nappi Code:1151640001

Elastic Adhesive Bandage
5 x 3M - EABBND5X3
Barcode :6009701456431
Nappi Code:1151641001

Elastic Adhesive Bandage
7.5 x 3M - EABBND75X3
Barcode :6009701456448
Nappi Code:1151642001

Elastic Adhesive Bandage
10 x 3M - EABBND10X3
Barcode :6009701456455
Nappi Code:1151643001

Elastic Adhesive Bandages

Safety Wear

MOBCP100W18




High Quality
Light Weight
Prevents Contamination
Double Stitched

Colour: White

Barcode:6009701456561
Nappi Code: 1155615001

Mob Cap White Double Stitched 100s

APRPE10015



Disposable Polyethylene Aprons

Pack Quantity:100s
Colour: White
- Oil, Liquid and Grease Resistant
- High Quality
- Economical and Disposable
- Tie Back

Pack: 100
White
Oil, Liquid and Grease Resistant
High Quality
Economical and Disposable
Tie Back

Barcode:6009701456615
Nappi Code:1155624001

**Apron 15 Micron -
121x72x0.013mm - 100s**

MOBCP100B18



High Quality
Light Weight
Prevents Contamination
Double Stitched

Colour: Blue

Barcode:6009701456578
Nappi Code:1155617001

Mob Cap Blue Double Stitched 100s

SHCVRNW100B



Overshoes specially adapted for health or hygienic purposes, e.g. comprising electrically conductive material allowing the discharge of electrostatic charges.

Disposable non-woven shoe covers pair Blue x 100 per pack

Barcode:6009701456608
Nappi Code:1155622001

Shoe Cover PE Blue 40x20cm - 100s

Safety Wear



Velcro Closure
Non – Woven
30gsm
Re-Usable

Small Barcode: 0710579687815 Medium Barcode: 0710579687808 Large Barcode: 0710579687822

Lab Coat - Non Woven White Velcro 30g

SLVCVR100W



Good Quality
Available in White
Packed in 100s
40x20cm

Barcode:6009701456585
Nappi Code:1155619001

**Sleeve Protector White PE 40x20cm
100s**

SLVCVR100B



Good Quality
Available in White
Packed in 100s
40x20cm

Barcode:6009701456592
Nappi Code:1155621001

**Sleeve Protector Blue
PE 40x20cm - 100s**

Wound Dressing



Gauze Swabs are used as a dressing in various applications primarily in wound care. Ideally designed to absorb and promote healing and preventing infection. Made from 100% high quality cotton exceeding industry.

Gauze Swabs 100 X 100mm 8 ply
GZCLIN100
Barcode :6009701452297
Nappi Code:1006029001

Gauze Swabs 75 X 75mm 8 ply
GZCLIN75
Barcode :6009701452303
Nappi Code:1006031001

Gauze Swabs 50 X 50mm 8 ply
GZCLIN50
Barcode :6009701453775
Nappi Code:1006033001

Gauze Swabs



- High quality
- Ventilated for quick healing
- Quality guaranteed FDA and CE approved
- Easy stretch and comfort
- Anti Allergy adhesive
- Protects and cushions wounds
- Sterile

Plasters Assorted 20s
PLAST20AS
Barcode :6009701452990
Nappi Code:1068534001

Plasters Fabric 20s
PLAST20FB
Barcode :6009701452976
Nappi Code:1068636001

Plasters Water Resistant 20s
PLAST20WR
Barcode :6009701452983
Nappi Code:1063841001

Plasters Fabric 100s pack
PLAST100FB
Barcode :6009701452846
Nappi Code:1068642001

Plasters Water Resistant 100s
PLAST100WP
Barcode :6009701452853
Nappi Code:1108750001

Plasters Kids Panda Forest 30s
PLAST30WR-K
Barcode :6009701451818
Nappi Code:1108811001

Plasters



Zinc Oxide Coating
Non Irritation Adhesive
Assists in faster healing
Gentle on Skin
Breathable and Flexible
Ideal for Dressing Fixation

50mm x 1m - FRP50MM
Barcode :6009701452945
Nappi Code:1184889001

25mm x 1m - FRP25MM
Barcode :6009701452938
Nappi Code:1184890001

12.5mm x 1m - FRP12MM
Barcode :6009701452921
Nappi Code:1184892001

Fabric Plaster Roll



A Skin like thin, sterile transparent film dressing that keeps out water, dirt, viruses and germs yet is breathable in the affected area. With a gentle adhesive on one side, it is impermeable to liquids and designed to be applied directly over clean skin or wounds.

- A waterproof, absorbent transparent dressing
- Sterile
- Latex free

- Hypoallergenic adhesive
- Breathable film
- All in one dressing, reducing application time

Dressing Transparent Film with
Frame style IV - 6x7cm U shape - 100s
TFDFSIV067
Barcode :6009701454000
Nappi Code:1135296001

Dressing Transparent Film with
Frame style IV - 7.2 x 9cm - 50s
TFDFSIV72X9
Barcode :6009701454017
Nappi Code:1135303001

Dressing Transparent Film with
Frame style IV - 9x11.6cm - 50s
TFDFSIV9116
Barcode :6009701453997
Nappi Code:1135304001

Dressing Transparent Film with Frame style IV

Wound Dressing



A Skin like thin, sterile transparent film dressing that keeps out water, dirt, viruses and germs yet is breathable in the affected area. With a gentle adhesive on one side, it is impermeable to liquids and designed to be applied directly over clean skin or wounds.

- A waterproof, absorbent transparent dressing
- Sterile
- Latex free
- Hypoallergenic adhesive
- Breathable film

Dressing Transparent Film Frame
Style 6 x 7cm - TDFD206X07 - 50s
Barcode :6009701453966
Nappi Code:1135305001

Dressing Transparent Film Frame
Style 10x 12cm - TDFD510X12 - 50s
Barcode :6009701453959
Nappi Code:1135306001

Dressing Transparent Film Frame Style



A Skin like thin, sterile transparent film dressing that keeps out water, dirt, viruses and germs yet is breathable in the affected area. With a gentle adhesive on one side, it is impermeable to liquids and designed to be applied directly over clean skin or wounds.

- A waterproof, absorbent transparent dressing
- Sterile
- Latex free
- Hypoallergenic adhesive
- Breathable film

Transparent Film with Non Adherent Pad -7x5cm
TFDWNAP7X5 - 50s
Barcode :6009701454093
Nappi Code:1135307001

Dressing Transparent Film with Non Adherent Pad -15x9cm
TFDWNAP15X9 - 50s
Barcode :6009701454116
Nappi Code:1135329001

Dressing Transparent Film with Non Adherent Pad -10x9cm
TFDWNAP10X9 - 50s
Barcode :6009701454109
Nappi Code:1135309001

Dressing Transparent Film with Non Adherent Pad - 25x10cm
TFDWNAP25X10 - 50s
Barcode :6009701454123
Nappi Code:1135330001

Dressing Transparent Film with Non Adherent Pad

PARGZE1010CH



Paraffin dressings provide a low-adherent first layer dressing that provides relief for the patient and facilitates the healing process.

Our open-weave gauze minimizes fraying and maintains shape.

Barcode :6009701455861
Nappi Code:1143427001

Paraffin Gauze 10cmx10cm - 10s



Non-woven Adhesive Wound Dressing is a sterile post-operative dressing indicated for the management of post-operative wounds and superficial wounds. It consists of soft non-woven fabric coated with one layer of hypoallergenic acrylic adhesive, an absorbent pad with low-adherent wound contact layer and two strips of release paper.

Water Resistant
Soft and comfortable adhesive
Latex Free

Non Woven with Non Adherent Pad - 6 x 8.3cm
NWDWNAP6X83 -50s
Barcode :6009701454048
Nappi Code:1135331001

Non Woven with Non Adherent Pad - 10 x 20cm
NWDWNAP1020 - 50s
Barcode :6009701454062
Nappi Code:1135333001

Non Woven with Non Adherent Pad - 10 x 30cm
NWDWNAP1030 - 50s
Barcode :6009701454079
Nappi Code:1135334001

Non Woven with Non Adherent Pad - 8 x 15cm
NWDWNAP8X15 - 50s
Barcode :6009701454055
Nappi Code:1135332001

Non Woven with Non Adherent Pad - 10 x 35cm
NWDWNAP1035 - 25s
Barcode :6009701454086
Nappi Code:1135335001

Dressing Non Woven with Non Adherent Pad

Tapes



Transparent, perforated film tape. Easy bi-directional tear. Allows the passage of water vapour from the skin. Provides excellent permeability to air and moisture.

Individually packed
Colour coded by size
Hypoallergenic
Dressing tape

Plastic Tape Transparent 1.25cm x 5m
PLTRNSP125

Barcode :6009701454413
Nappi Code:1135271001

Plastic Tape Transparent 2.5cm x 5m
PLTRNSP250

Barcode :6009701454437
Nappi Code:1135282001

Plastic Tape Transparent 5cm x 5m
PLTRNSP500

Barcode :6009701454451
Nappi Code:1135283001


Plastic Tape Transparent 10cm x 5m
PLTRNS1000

Barcode :6009701454499
Nappi Code:1135293001

Plastic Tape Transparent 7.5cm x 5m
PLTRNSP750

Barcode :6009701454475
Nappi Code:1135287001

Plastic Tape Transparent



Chemically non reactive and biologically inert
Non-woven Paper
Latex-Free
Hypoallergenic
Excellent air permeability
Extremely easy tear by hand
Available in various sizes

Microporous Paper Tape 1.25cm x 5m
PMCRPT125

Barcode :6009701456462
Nappi Code:1151644001

Microporous Paper Tape 5cm x 5m
PMCRPT500

Barcode :6009701456486
Nappi Code:1151647001

Microporous Paper Tape 10cm x 5m
PMCRPT1000

Barcode :6009701456547
Nappi Code:1151650001

Microporous Paper Tape 2.5cm x 5m
PMCRPT250

Barcode :6009701456479
Nappi Code:1151645001

Microporous Paper Tape 7.5cm x 5m
PMCRPT750

Barcode :6009701456493
Nappi Code:1151649001

Microporous Paper Tape



Zinc oxide plaster tape is a non-stretch adhesive tape. Its usages are various; they are often used to prevent injuries, protect wounds or heal cuts faster.

Quality non stretch plaster

Extra Strong adhesive

Hospital pack with rolls of 6 or 12 per box depending on size

Zinc Oxide Plaster Tape 1.25x5m
ZINCOXD125

Barcode :6009701454673
Nappi Code:1135253001

Zinc Oxide Plaster Tape 2.5x5m
ZINCOXD250

Barcode :6009701454697
Nappi Code:1135254001

Zinc Oxide Plaster Tape 5x5m
ZINCOXD500

Barcode :6009701454680
Nappi Code:1135255001

Zinc Oxide Plaster Tape 10x5m
ZINCOXD1000

Barcode :6009701454710
Nappi Code:1135257001

Zinc Oxide Plaster Tape 7.5x5m
ZINCOXD750

Barcode :6009701454703
Nappi Code:1135256001

Zinc Oxide Tape

Air Management



The Guedel Airway is a Oropharyngeal Tube used to prop open the upper airways. Due to it's form, the correctly used Guedel Airway prevents that the tongue of a patient can drop back and therefore block the airway. It does not offer protection from Aspiration. The use of a suction catheter is however possible. Furthermore, the Guedel Airway can be used for tracheal tube fixation. For this use, the Guedel Airway also functions as bite block to prevent an occlusion of the tracheal tube caused by the patient's teeth

Guedel Airway Number 00 - 50mm Blue
GUEAIR00
Barcode : 0710579687600

Guedel Airway Number 0 - 60mm Black
GUEAIRW0
Barcode : 0710579687617

Guedel Airway Number 1 - 70mm White
GUEAIR01
Barcode : 0710579687624

Guedel Airway Number 2 - 80mm Green
GUEAIR02
Barcode : 0710579687631

Guedel Airway Number 3 - 90mm Yellow
GUEAIR03
Barcode : 0710579687648

Guedel Airway Number 4 100mm - Red
GUEAIR04
Barcode : 0710579687655

Guedel Airway Number 5 110mm
- Light Blue
GUEAIR05
Barcode : 0710579687662

Spirometer



CHOXCN

Continuous Mode 1 - 5L Adjustable
Variable frequency compressor
Vehicle Power Supply

Barcode: 0710579687570
Nappi Code: N/A

Oxygen Concentrator with Battery



Neck Brace

Adult: CHNBAD
Barcode : 0710579687587
Nappi Code:

Child: CHNBCH
Barcode : 0710579687594
Nappi Code:

Neck Brace

Air Management



Nebulization kit suitable for use with most nebulizer models
Includes one face mask, tubing and a medication cup

Adult: NEBMACH Child: NEBMCCH
Barcode:6009701459944 Barcode:6009701459951
Nappi Code:1184896001 Nappi Code:1184897001

Nebulizer Mask



Easily transportable from one hospital ward to another, particularly suitable for tracheotomized patients, minor surgical applications and postoperative therapy at home
To be used for aspirating body fluids such as mucus.

Barcode: 6009701459722

Sputum Aspirator - MY696



Oxygen Rebreather Mask
Packed in 100s

Adult: CBMSKRBA
Barcode: 0710579687792
Nappi Code: N/A

Oxygen Rebreather Mask - Adult



A Spirometer is helpful in assessing conditions such as asthma, pulmonary fibrosis, cystic fibrosis and COPD. Just inhale through Lung Exerciser for a few minutes twice a day in a simple breathing pattern as shown in the instructions.

Barcode: 6009701459524

Spirometer

Nasal Cannula with 1.5m supply tube
Single Patient use.
Do not sterilize
Over the ears

Nasal Cannula - Adult
NASCANAD1CH
Barcode :6009701455830
Nappi Code:114348001

Nasal Cannula - Child
NASCANCHD1CH
Barcode :6009701455847
Nappi Code:114340001

Nasal Cannula - Neonatal
NASCANNEO1CH
Barcode :6009701455847
Nappi Code:114340001



Nasal Cannula

Needles and Syringes



High Quality
Sterile
Non Pyrogenic
Latex Free
Manufactured under ISO
supervision
CE Approved

Syringe 3ML + 22G Needle-100s (L/S)
SYR3ML22G
Barcode :6009701452822
Nappi Code:1129538001

Syringe 3ML + 21G Needle-100s (L/S)
SYR3ML21G
Barcode :6009701451696
Nappi Code:1129539001

Syringe 5ML + 22G Needle-100s (L/S)
SYR5ML22G
Barcode :6009701452839
Nappi Code:1047741001

Syringe 1ML (L/S) + 27G 1/2" (40X13)
SYR1ML27G
Barcode :6009701453195
Nappi Code:1087696001

Needles and Syringes - Combo Packs



Single Use Only
Dispose used syringes into bio-hazard containers
For use by trained personnel
Follow infection control procedure during preparation & administration

SYR1MLLS
Syringe 1ml (L/S) - 100s

Barcode:6009701450040
Nappi Code:1092561001

SYR2MLLS
Syringe 2ML (L/S) - 100s

Barcode:6009701450866
Nappi Code:1087676001

SYR3MLLS
Syringe 3ML (L/S) - 100s

Barcode:6009701450279
Nappi Code:1068190002

SYR5MLLS
Syringe 5ML (L/S) - 100s

Barcode:6009701450286
Nappi Code:1068191001

SYR10MLLS
Syringe 10ML (L/S) - 100s

Barcode:6009701450293
Nappi Code:1068195001

SYR20MLLS
Syringe 20ML (L/S) -50s

Barcode:6009701450541
Nappi Code:1068196001

SYR50MLLS
Syringe 50ml 3 Part L/Slip 25s

Barcode:6009701457247
Nappi Code:1155943001

Syringes

Needles and Syringes



High Quality
Sterile
Manufactured in an ISO approved facility,
CE Approved
Packed in 100's

NDL18GCH
Needle 18G 1&1/2" (1.2x38) - Pink

Barcode:6009701450668
Nappi Code:1068644001

NDL21GCH
Needle 21G 1&1/2" (80x38) - Green

Barcode:6009701452396
Nappi Code:1005667001

NDL21G8025CH
Needle 21G 1" (80x25) - Green

Barcode:6009701453218
Nappi Code:1087692001

NDL22GCH
Needle 22G 1 & 1/4" (.70x32) - Black

Barcode:6009701452402
Nappi Code:1005669001

NDL22G7025CH
Needle 22G 1" (.70x,25) - Black

Barcode:6009701453225
Nappi Code:1087688001

NDL23GCH
Needle 23G 1 & 1/2" (64x38) - Blue

Barcode:6009701450255
Nappi Code:1068192001

NDL23G6025CH
Needle 23G 1" (.60X,25) - Blue

Barcode:6009701453232
Nappi Code:1087680001

NDL25GCH5016
Needle 25G X 5/8" (50X16) - Orange

Barcode:6009701450606
Nappi Code:1092562001

NDL25GCH1,0
Needle 25g 1 inch (0.50x25mm) Orange

Barcode:6009701451344
Nappi Code:1096991001

NDL25GCH1,5
Needle 25g 1 1/2 inch (0.50 x38mm) Orange

Barcode:6009701451375
Nappi Code:1096993001

NDL26G4513CH
Needle 26G 1/2" (0.45X13) - Brown

Barcode:6009701453201
Nappi Code:1087684001

NDL27G4013CH
Needle 27G 1/2" (40X13) - Grey

Barcode:6009701453249
Nappi Code:1087686001

NDL27G4013CH
Needle 27G 1/2" (40X13) - Grey

Barcode:6009701453249
Nappi Code:1087686001

INS29G05ML
Insulin Needle 29G x 0,5ML

Barcode:6009701450637
Nappi Code:1087714001

INS29G1ML
Insulin Needle 29G x 1ML

Barcode:6009701450644
Nappi Code:1068193001

SCALVEIN21GCH
Scalp Vein Set Continuous - 21G

Barcode:6009701458770
Nappi Code:1166334001

SCALVEIN23GCH
Scalp Vein Set Continuous - 23G

Barcode:6009701450682
Nappi Code:N/A

NDLAC3025CH
Needle Acupuncture .30x25mm

Barcode:60097014585541
Nappi Code:1166267001


NDLAC3525CH
Needle Acupuncture .35x25mm

Barcode:6009701458558
Nappi Code:1166267001

Needles

Needles and Syringes

SPN1890CH




The defining feature of a Quincke needle is its sharp, beveled tip.

This design allows the needle to cut through the dural and arachnoid membranes with greater ease than needles with a rounded, non-cutting tip.

Barcode: 0710579688300

Spinal Needle Quincke Point 18gx90cm - 50s

SPN2090CH




The defining feature of a Quincke needle is its sharp, beveled tip.

This design allows the needle to cut through the dural and arachnoid membranes with greater ease than needles with a rounded, non-cutting tip.

Barcode: 0710579688317

Spinal Needle Quincke Point 20gx90cm - 50s

SPN2290CH




The defining feature of a Quincke needle is its sharp, beveled tip.

This design allows the needle to cut through the dural and arachnoid membranes with greater ease than needles with a rounded, non-cutting tip.

Barcode: 0710579688324

Spinal Needle Quincke Point 22gx90cm - 50s

SPN2390CH




The defining feature of a Quincke needle is its sharp, beveled tip.

This design allows the needle to cut through the dural and arachnoid membranes with greater ease than needles with a rounded, non-cutting tip.

Barcode: 0710579688331

Spinal Needle Quincke Point 23gx90cm - 50s

SPN2590CH



The defining feature of a Quincke needle is its sharp, beveled tip.

This design allows the needle to cut through the dural and arachnoid membranes with greater ease than needles with a rounded, non-cutting tip.

Barcode: 0710579688348

Spinal Needle Quincke Point 25gx90cm - 50s



CLINIhealth
leaders in healthcare

Gloves



High quality Nitrile gloves -Latex free offers superior protection
ISO certified
Use in General Medical and dental industries
Packed in 100s / weight

Nitrile gloves provide superior quality

- Allergy free
- Latex free
- ideal for home and professional
- Powder free
- Product packaging may differ

Nitrile Gloves - Small-100s
GLNCHS

Barcode :6009701452709

Nappi Code:1028258001

Nitrile Gloves - Medium-100s
GLNCHM

Barcode :6009701452716

Nappi Code:1028260001

Nitrile Gloves - Large-100s
GLNCHL

Barcode :6009701452723

Nappi Code:1028262001

Nitrile Gloves Blue



High Sensitivity and durability
Works well in high risk operations
Cost effective and chemically resistant
Powder free
Biodegradable
Textured Grip
Made with high grade rubber

Latex Gloves- Small 100s
GLPFCHS

Barcode :6009701452679

Nappi Code:1027489001

Latex Gloves- Medium 100s
GLPFCHM

Barcode: 6009701452662

Nappi Code:1027487001

Latex Gloves- Large 100s
GLPFCHL

Barcode :6009701452655

Nappi Code:1027488001

Latex Gloves - White

High quality Nitrile gloves (Latex free) offers superior protection
ISO certified
Use in General Medical and dental industries
Packed in 100s / weight
Allergy free
Latex free
Ideal for home and professional
Powder free
Product packaging may differ

Nitrile Gloves Black
Small 100s
GLNBCHS

Barcode :6009701451412

Nappi Code:1119321001

Nitrile Gloves Black
Medium 100s
GLNBCHM

Barcode :6009701451429

Nappi Code:1119323001

Nitrile Gloves Black
Large - 100s
GLNBCHL

Barcode :6009701451436

Nappi Code:1119326001



Nitrile Gloves Black

Gloves

PE (Polyethylene)
Multi Purpose Gloves
Measurement: One Size Fits
Colour: Clear
100 in a pack

GLPLST100LCH



Barcode :6009701455878

Nappi Code:N/A

Disposable PE Gloves - 100s



High quality Nitrile gloves -Latex free offers superior protection
ISO certified
Use in General Medical and dental industries
Packed in 100s / weight

Nitrile gloves provide superior quality

- Allergy free
- Latex free
- ideal for home and professional
- Powder free
- Product packaging may differ

Nitrile Gloves Pink - Small-100s
GLNPCHS

Barcode :710579687105

Nappi Code:1216945001

Nitrile Gloves - Medium Pink-100s
GLNPCHM

Barcode :710579687112

Nappi Code:1216951001

Nitrile Gloves Pink

PE (Polyethylene) – Multi Purpose Gloves
Measurement

One Size Fits
Colour

Clear
50 in a pack

CHEBPES0



Barcode : 0710579687785

Nappi Code:N/A

PE Gloves – Shoulder Length 50s

First Aid Kits

FIRAICOM



First aid protection – To keep you safe from minor injuring when travelling
Contained in a compact bag – Durable, value for money and easy to pack
Size: 110mm x 75mm x 45mm for your convenience


Barcode :N/A
Nappi Code:N/A

First Aid Kit Compact

34 Piece First Aid Kit Contents

- 1 Pair (2) x Examination Gloves Large
- 1 x Safety Scissors
- 2 x Cleansing Alcohol Wipes
- 3 x Insect Sting Relief Pad
- 4 x Alcohol Prep Pads
- 2 x Antiseptic Wipes
- 5 x Plasters 76x45mm
- 10 x Plasters 76x19mm
- 2 x Non Woven Gauze Swabs 5x5cm
- 2 x Sterile Gauze 5x5cm
- 1 x Bag

FIRB34PCH



Barcode :6009701450163
Nappi Code:N/A

First Aid Bag 34Pc

Regulation 3
First Aid Kit (Government Specs)
High quality
ISO certified
Ideal for factory offices and homes
In Red Bag

FAKBG30CH



Barcode :6009701458862
Nappi Code:N/A

First Aid Kit Regulation 3 in Red Bag

Motorist
Contents List

- 1 x A4 Bag
- 1 x Eye Pad
- 1 x Triangular Bandages – non woven
- 2 x First Aid Dressing #3
- 1 x First Aid Dressing # 5
- 2 x Conforming Bandage 75mm
- 1 x Rescue Blanket

FAKBGMTCH




Barcode :6009701458473
Nappi Code:N/A

First Aid Kit Motorist in Red Bag

Government Regulation 7
First Aid Kit

FAK701PK




Barcode :N/A
Nappi Code:N/A

First Aid Kit Regulation 7 in Metal Box

Government Regulation 3
First Aid Kit

FAK301PK



Barcode :N/A
Nappi Code:N/A

First Aid Kit Regulation 3 in Metal Box

First Aid Kits

FACLI310



First Aid Kit with 303 essential first aid items.

Comes in zipper bag

Email us for full list of contents

Barcode :6009701450767
Nappi Code:N/A

First Aid Kit 303 Pieces in Zip Bag Clinihealth

FACLI210




First aid kit containing 210 items

In convenient zipper bag

Antiseptics
Pack of 100 strips

Barcode :6009701450767
Nappi Code:N/A


First Aid Kit 210 Pieces in Zip Bag CliniHealth



PVC	SILICONE
AMBUPVACH - Adult Barcode :6009701459418	AMBUPVICHs - Infant Barcode :0710579687396
AMBUPVCCH - Child Barcode :6009701459425	AMBUPVCCHs - Child Barcode : 0710579687402
AMBUPVICH - Peads Barcode :6009701459432	AMBUPVACHs - Adult Barcode : 0710579687419

Oxygen Resuscitator Ambu Bag

CHCPRBR




High density polyethylene plastic
Size: 49x41x8.5cm

Barcode: 0710579687679

CPR Board

CLIEAR100



Head Immobilizer

Barcode :
Nappi Code:N/A

Head Immobilizer

Cotton Products

CLIEAR100




High quality Pure Cotton earbuds
Soft and gentle on skin
Use for medical and personal skincare routine
Safe and gentle
Better pricing on bulk orders
100 earbuds

Barcode :6009701450477
Nappi Code:N/A

Earbuds Luxury 100 Clinihealth

CLICOT25G



High quality Pure Cotton Roll
Soft and gentle on skin
Use for medical and personal skincare routine
Safe and gentle
Better pricing on bulk orders

Barcode :6009701450484

Cotton Roll Luxury Soft 100g CLICOT100G
Barcode :6009701450514

Cotton Roll Luxury Soft 500g CLICOT500G
Barcode :6009701450521

Cotton Rolls

COTT500G



High quality Pure Cotton ball
Soft and gentle on skin
Use for medical and personal skincare routine
Safe and gentle

Barcode :6009701452365
Nappi Code:N/A

100G Balls - COTT100G
Barcode :6009701452761
Nappi Code:N/A

Cotton Balls

CLICOTR80



100% pure, soft and gentle cotton
Dual textured
Ideal for application of lotions, creams, and make-up
Reverse side provides resistance required for removal of makeup and cleaning of the skin
Easy to use drawstring closure bag for convenient and re-use

Barcode :6009701451610
Nappi Code:N/A

Cotton Rounds 80s

EYEPDSS50



Non adhesive absorbent eye pad is made of an elastic non-woven dressing with absorbent pad. 8 x 6cm
Individually packed in sterile, peel open pouches
Designed to protect wound site while promoting blister pouch
Non-stick outer facing allows the dressing to absorb, but not interfere with wound

Barcode :6009701453980
Nappi Code:1135276001

Eye Pad Sterile with non adherant pad - 50

Linen Savers



Breathable Top Sheet
Special Embossed Surface
Increased Absorbing Effect
Super Absorbent – Turns liquid into Gel
Soft Waterproof Backsheet
Blocks Unpleasant odour

Linen Savers 60x60 -Retail 20s
LIN208PL6060

Barcode :6009701451931

Nappi Code:1087744001

LINEN SAVERS 60x60 (200s)
LIN200SPL6060

Barcode :6009701451955

Nappi Code:1087744002

Linen Savers 60x60



Breathable Top Sheet
Special Embossed Surface
Increased Absorbing Effect
Super Absorbent – Turns liquid into Gel
Soft Waterproof Backsheet
Blocks Unpleasant odour

LINEN SAVERS 60x90 -RETAIL 20s
LIN208PL6090

Barcode :6009701451948

Nappi Code:1087761001

LINEN SAVERS 60x90 (200s)
LIN200SPL6090

Barcode :6009701451962

Nappi Code:1087761002

Linen Savers 60x90

Adult Diapers



Product Details
Instant Absorption
Consist of mixed polymer absorbent and villous pulp which results to fast absorption, instant drying, leakage prevention, restraint odour, reduced bacterial growth, and provide the feeling of freshness.

M - 71x120cm
L - 96x140cm
XL - 111x160cm

Adult Diapers Premium - Medium
ADUDIAPM10

Barcode :6009701455526

Nappi Code:1135277001

Adult Diapers Premium - Large
ADUDIAPL10

Barcode :6009701454994

Nappi Code:1135278001

Adult Diapers Premium - XL
ADUDIAPXL10

Barcode :6009701455007

Nappi Code:1135279001

Adult Diapers Premium - Pull Ups



Great absorbent, vacuum package
Super absorbent
Tear resistant top sheet.
Scientifically developed, super absorbent core to quickly lock in wetness.
115x220cm Mattress

Barcode :710579687532
Nappi Code:N/A

Bed Sheet Non Woven 115x220mm



Product Details
Instant Absorption

Consist of mixed polymer absorbent and villous pulp which results to fast absorption, instant drying, leakage prevention, restraint odour, reduced bacterial growth, and provide the feeling of freshness.

Adult Diapers Premium Strap On
Medium - ADUDIASTM14

Barcode :6009701457278

Nappi Code:1149550001

Adult Diapers Premium Strap On
Large - ADUDIASTL12

Barcode :6009701457285

Nappi Code:1149551001

Adult Diapers Premium Strap On
Xlarge - ADUDIASTXL10

Barcode :6009701457292

Nappi Code:1149552001

Adult Diapers Premium - Strap Ons

Measuring and Checking Apparatus



Easy to Read: High Quality quartz movement nurse, Unique safety Clip-on pin design fob watch.

This model features a clear white dial with black Arabic numerals and red pulsation numbers around the edge.

Nurses Watches Stainless Steel
WATNURSST

Barcode :6009701450859

Watch Nurses (white Silica Gel)
WATNURWSG

Barcode :6009701450842

Watches



MEACUCCH

- Premium 250ml glass conical measuring cup
- Ungraduated
- Pouring spout makes for convenient and accurate distribution of liquids
- Comes with a stable yet elegant rounded base

Barcode 6009701458879

Measuring Cup Conical 250ml



SPECT20CH

Easy to read graduations on two sides in 10ml increments
Red polyethylene screw cap with ribbed edge for extra grip for effective sealing and easy opening
Expiry date, Lot Number and sterile warning are clearly printed on each wrapped container

Barcode :6009701458886

Specimen Container Sterile 120ml



PILCNTCH1

A pill counting tray is a quick and sanitary way to handle and count pills. Pharmacy and home users can use it to verify prescription count accurately, avoid miscounted divide plenty of pills tablet, even capsule, into smaller increments. Also can count/sort small items

Barcode :6009701452624

Pill Counter with Spatula



ANESPHYCH

- Markings on cuff to ensure proper positioning
- Cuff includes convenient gauge holder
- Standard cotton cuff with Velcro closure
- No-pin stop manometer (Gauge) 300mmHg

Barcode : 6009701450989

Aneroid Sphygmometer



ELHEIWEICH

Max Capacity : 250kg
Stature Meter: 70-190cm
Metal Height Measure Rod
Max Weight: 200kg
Min Weight Division: 0.1kg
Height Range: 70-190cm
Min Height Per Division: 5mm
Load Platform Area (L x W): 375x275mm

Barcode :6009701457964

Nappi Code:N/A

Electronic Height and Weight Scale



ELBBYSCCH

This balance is specially made for measuring the weight of a baby
The pan may be removed to weigh babies under 20kg standing up
This instrument is ideal for a family, hospitals or any infant care center
Kg/lb switchable

Barcode :6009701456974

Electronic Baby Scale



Arm Circumference MUAC Measuring Tape for children and adult, graduation in centimeter and capable to use for taking of measurement of MUAC.
Weather-resistant, flexible, non-tear synthetic single-slotted insertion tape; approved color-coded division for mid-upper arm circumference 26 cm capacity

Adult MUAC Tapes are also standardized by WHO and UNICEF with colour codes.
Breathable film

Muac Tape - Adult 56.5 cm

MUACADCH

Barcode :6009701456981

Muac Tape - Child 26.5cm

MUACCHCH

Barcode : 6009701456998

Muac Tape

Single Use Packed in convenient resealable aluminium pouch Latex and PVC free ECG Electrodes with pregelled AG / AGCL sensor



ECGPRADCH

Barcode :6009701457957

ECG Electrodes - Pre-Gelled Adult

Measuring and Checking Apparatus

BBSCDIYCH



Weighs up to 20 kg
No batteries required.
Plastic bucket seat for infant to lie down.

Barcode :6009701457018
Nappi Code:N/A

Baby Scale with Dial - Yellow 20kg

BBSCDIBWCH



Weighs up to 20 kg
No batteries required.
Plastic bucket seat for infant to lie down.

Barcode: 6009701457025
Nappi Code:N/A

Baby Scale with Dial - Blue and White 20kg

CHBSM



Bathroom Scale(mechanical)
Capacit:140kgs(300lb)
Grad:1kg(1lb)
Dimension:28.7x26.7x5.5cm

Barcode: 600970145944

Bathroom Scale - Mechanical



Small 22-32cm
Medium - 22 - 36cm
Large - 22 - 42cm

Barcode

Small: 6009701459456
Medium: 6009701459463
Large: 6009701459470

CHBPCFS / CHBPCFM / CHBPCFL

BP Cuff - Small / Medium / Large

CHSGTIM




Designed with enough large LCD display, you can easily read the digits, very convenient
Loud alarm sound, you can hear the alarm in another room
With a magnet in the clip of the timer, you can attach timer to any item made of iron and steel.

Barcode : 6009701459494

Single Group Timer

CHPK01B



Instant readings help you track any changes in your airflow
Compact design allows you to check your PEFR discreetly any time you want
Assists detection of a pending asthma attack

Barcode :6009701459517

Peak Flow Meter Blue

CH4GPTIM



Built in stand for use on the counter
Magnet for easy attachment to coffee machines and other appliances
Batteries included
4 separate timers / groups
Ability to count up or count down to an alarm
Easily set hours, minutes and seconds separately

Barcode :6009701459500

4 Group Timer

Measuring and Checking Apparatus



DGS200LP

Weight Range : 1-200kg
Function : Body Weight, BMI
Display Type : LCD
Power: Rechargeable Battery
Shape : Polygon
Suitable for: Body Check
Certification : CE, RoHS, ISO
Appearance Screen Color
Color: Black
Capacity : 200kg

Barcode :0710579687310

Digital BMI Scale TCD200LP



DCPUPTCH

Material: stainless steel
Light source: LED light source
Power supply: rechargeable lithium battery
Battery capacity: 70 mA
Color rendering index: greater than 90

Barcode :6009701456950

Doctors Pupil Torch Rechargeable

A stethoscope helps you to listen to the internal sound of the heart, lungs and bowels.



STETWTIMCH

Barcode :6009701457957

Stethoscope With Timer



CHSP5025

Screw on green cap Urine Cup, x 25

Sampling bottle.

A cylindrical container, usually closed at a chosen depth, to trap a water sample and transport it to the surface without introducing contamination.

Barcode: 6009701459586

Specimen Jars 50ml - 25s



MCHEIWEICH

Height and Weight scale
All types of height and weight mechanical and electronic scales importers and dealers in the market.
Max weighing: 180kg
Range to be measured: 850-2100mm
Min value of height per division: 1mm
W.: 15kg N. W.: 13kg

Barcode :6009701459760

Nappi Code:N/A

Mechanical Height and Weight Scale



CHCS250

Dimensions: 130x69x93cm
With rechargeable battery (included)
Foldable handle
Large LCD displayer 5.2inch with blue backlight
360degree turning wheels
Net weight:23.2kg
Gross Weight: 34.2kg

Barcode :6009701459753

Nappi Code:N/A

Chair Scale - Load Capacity 250Kg

Measuring and Checking Apparatus



CHLWM2
Height Measure 2m Wall Mounted
Sold Individually
Lightweight and Portable
Barcode : 0710579687037
Nappi Code:N/A

Wall Mounted Height Measure




CHPTSTA
A stadiometer is a piece of medical equipment used for measuring human height.
It is usually constructed out of a ruler and a sliding horizontal headpiece which is adjusted to rest on the top of the head.
Barcode: 0710579687051

Portable Height Stadiometer - 20 - 210cm



CHMNM102
Large LCD display
90 groups of memory
Indication for irregular heartbeats
Detection for cuff wearing
Indication for Movement Error
Average of three measurements
Rechargeable lithium battery
Indication for movement error
Detection for cuff wearing
Barcode : 6933257922365

Mesh Nebulizer - Yuwell M102



CHPK01B
Instant readings help you track any changes in your airflow
Compact design allows you to check your PEFR discreetly any time you want
Assists detection of a pending asthma attack
Barcode : 6933257915510


Basic Blood Pressure Yuwell YE8600D



CH4GPTIM
Built in stand for use on the counter
Magnet for easy attachment to coffee machines and other appliances
Batteries included
4 separate timers / groups
Ability to count up or count down to an alarm
Easily set hours, minutes and seconds separately
Barcode : 6933257913004

Basic Blood Pressure Yuwell YE-660D

Measuring and Checking Apparatus



CHTWTH

Specification
Temperature range
Indoor: -10 degrees Celsius to 50 degrees Celsius
Outdoor -50 to +70 degrees Celsius Temperature range
Indoor: -10 degrees Celsius to

Barcode: 6009701459913

Twin Reading Fridge Thermometer -with Probe




CHIETHRYHT101

Measurement range:32-43
Temperature unit
Dimension:149*38*42mmResp
onse time: 1 second Auto shut
off: 60 seconds Measure
distance: 0-5cmMeasuring
Color Indicator

Barcode: 6933257922082

Infrared Ear Thermometer Yuwell- YHT101



CHIETHRYHT101

Patella hammer
St Steel
Babinsky

Barcode: 710579687280

Babinsky Patella Hammer



CHAS160

Max weighing capacity: 150kg.
Black rubber top.

Barcode: 710579687044

Adult Scale Max 160kg




MEATUM30

Medicine Tumbler 50ml Plastic
-50s

Barcode: 0710579687426

Medicine Tumbler 30ML




CHEACKLT

Ear Cleaning Kit

Barcode: 0710579687686

Ear Cleaning Kit




ANESPHSTCH

Markings on cuff to ensure proper positioning
Cuff includes convenient gauge holder
Standard cotton cuff with Velcro closure
Comes with zipper case
Includes Stethoscope

Barcode: 710579687297

Aneroid Sphygmomanometer with Stethoscope



CHLMEPK

Ear Cleaning Kit

Barcode: 0710579687693

Luminous Ear Pick

Measuring and Checking Apparatus

CHBMIC




The Body Mass Index (BMI) is one way to measure body size. It's a tool to estimate body fat and screen for obesity and health risks. It can be calculated with a BMI calculator and classifies people as being underweight, overweight, and obese based on their height and weight

Barcode :6009701459531

BMI Calculator 1.5M

CHPRWH




PVC Pregnancy Wheel
Diameter:108mm

Barcode: 6009701459548

PVC Pregnancy Wheel

CHRTHE



Barcode :6009701459487

Room Thermometer / Hygrometer

FIRDEFRCH



Feature:

- Two-button operation
- Extensive voice and visual prompts for the operator
- Biphasic energy output
- Lock-out protection to prevent inadvertent defibrillation
- Device has an internal memory(60 hours)
- ECG signals and events can be transferred to the computer.
- It has a Long-life disposable LiMnO2 type battery (at least 5

Barcode :6009701459555

5MLDOSCH




Dosing syringe 0.5 ml – 5 ml.
Smallest unit of measurement 0.5 ml.

Barcode: 0710579688560

5ml Dosage Syringe Non Sterile - Individually Packaged

MEADP5CH




Medicine Dropper

Barcode: 0710579688577

Medicine Dropper 5ml Plastic

FTSTCPCH



A device used to obtain information about a fetus within the uterus.

Aluminum pinard is made of high-quality material.

Barcode:0710579688584

Fetal Stethoscope

CHDEF5S



Biphasic energy delivery
Manually switchable language

Barcode :0710579688591

AED First Aid Defibrillator Defi 5S

Measuring and Checking Apparatus

CHDRBT15/30/50



Amber brown glass bottle with a dropper lid
Perfect for storing essential oils, perfumes, pharmaceuticals and many other liquids
Ideal for storing, decorating and organizing.
Made of durable glass construction

15ml Barcode : 0710579688607
30ml Barcode : 0710579688614
50ml Barcode : 0710579688621

Dropper Bottles Amber Glass 15/30/50ml

CHOIJA200/100/50



PVC Pregnancy Wheel
Diameter:108mm

200g Barcode: 0710579688638
100g Barcode: 0710579688645
50g Barcode: 0710579688652

Ointment Jar Glass Amber

CHTABVP20/50/75/100



Tablet Vials Plastic Snap on Lid Clear

20ml Barcode : 0710579688669
50ml Barcode :0710579688676
75ml Barcode :
100ml Barcode :

Tablet Vials Plastic Snap on Lid Clear

Packaging and Storage

TABBAG1000



Re-sealable clear plastic bag, made of Pure Low Density Polyethylene (LDPE).
Bag Dimensions: 115mm x 75mm
Zip lock.
Easy to use white bag with name expiry and batch number
Packed in 1000s
Certified for medical use

Barcode :6009701451122

Tablet Bags - 1000s

PILLBOXTDCH1



The Seven-Day Pill Organizer is perfect for helping you organize your pills.


Seven individual compartments allow you to fill in a week's medication at one time.

Each compartment can hold small, daily doses of medication or vitamins.

Barcode :6009701455823
Nappi Code:1143431001

7 Day Pill Box

PILLGR14CH



Morning and Noon Pill organizer
With 14 Grids

Barcode :6009701456929
Nappi Code:1155947001

Pill Organizer 14 Grid - Morn & Noon



MOBILITY



Bathroom Safety



CH7952

Multi-Functional Bathroom Bench shower stool can be installed at bathroom, stairs, kids room or any room in the house. Sturdy Construction

Folding shower seat - heavy duty

- Seat Width: 35cm
- Seat Height : 36cm
- Maximum Weight Capacity: 130kg

Barcode:6009701456073

Nappi Code:1143636001

Shower Chair - Wall Mounted Foldable



CH568

Single Stepper for Bath or Shower makes getting in and out of the bath tub easy and safe.

It is also ideal for those who need a bit of help getting in and out of the bed.

A versatile and lightweight item.

- Width: 41cm
- Height: 15cm
- Length: 25cm
- Net Weight: 3.2kg

Barcode:6009701453348

Nappi Code:1103973001

Single Stepper Bath/Shower



CH569

Double Stepper for Bath or Shower makes getting in and out of the bath tub easy and safe.


A versatile and lightweight item.

- Total Width: 41cm
- Max Weight: 100kg
- Gross Weight: 3.36kg
- Nett Weight: 2.64kg

Barcode:6009701451856

Nappi Code:1103401001

Double Stepper Bath/Shower



Bath Bench with Back rest. Comfortable and convenient lightweight Bench that is easy to apply to bath and remove.

- Seat Width: 41cm
- Total height: 46cm
- Length: 39cm
- Max Weight: 120kg
- Gross Weight: 3.2kg
- Nett Weight: 2.5kg

Barcode:6009701451894

Nappi Code:1103405001

Bath Bench with Backrest



CH794LBS



CH7946

L-Shaped Grab Rail for the shower can be positioned on any flat wall and comes in stainless steel, it is designed for elderly and disabled persons

The L-shaped / Shower Grab Rail provides a sturdy hand hold with anti slip rubber grip.

Net Weight: 5.82kg

Barcode:6009701453379

Nappi Code:1103903001

Hand Rail L Shaped

Grab rails next to the toilet, bath or basin will ensure you are safe in these areas.

- Small: 30cm
- Medium: 40cm
- Large: 49cm
- XL: 50cm
- XXL: 60cm

Barcode:6009701453386

Nappi Code:1103458001

Grab Rail - Small / Medium / Large / XL / XXL



CH7947S
CH7947M
CH7947L
CH7947XL
CH7947XXL

Bathroom Safety



CH7948

Support yourself easily and safely with the Drop Down Rail With Leg. This rail can fit onto the toilet wall and features a drop-down leg for additional support.

- Net Weight: 2.35kg
- Gross Weight: 2.44kg

Barcode:6009701453379
Nappi Code:1103903001

Drop Down Rail with Leg



CH7944L

Handicap grab bars that allow you to sit safely or stand at the toilet. The toilet stand is easy to install and setup. The toilet rail would pair up nice with a handicap toilet seat for added bathroom safety.

- Width: 41cm
- Height: 15cm
- Length: 25cm
- Net Weight: 6.5kg

Barcode:6009701451979
Nappi Code:1103462001

Toilet Safety Rails




CH671B

N.W.:1.7kg G.W.:2.2k
Description
Raises the seating height from 5 to 15 cm
Height adjustable
Firmly secures onto toilet bowl
Suitable for hip replacement patients
Weight capacity 100kgg

Barcode: 6009701459265
Nappi Code:

Toilet Seat Raiser without handles




CH670B

Raised toilet seat. Easily attaches to a toilet. Aids user in sitting and rising. Circumference - 25 - 37.3cm Length - 20cm

Barcode:6009701451979
Nappi Code:1103462001

Toilet Seat Raiser with Handles



CHX00201

Foldable shower seat wall mounted, saves space
Steel frame durable, rubber tip, height adjustable.

Seat Width: 51.5 cm
Upholstery: PE
Weight Capacity: 130 kg
Nett Weight: 4.75 kg
Overall Width: 51.5 cm

Barcode:6009701456172
Nappi Code:1143637001

Shower Chair- Wall mounted with legs



CH797LH

This bath chair shower seat assembles in just minutes with no tools required.

- Total Height: 39-49cm
- Total length: 39cm
- Seat Width: 50.5cm
- Max weight: 160kg
- Length Total:46cm
- Seat Height: 36 – 46cm

Barcode:6009701455540
Nappi Code:1143638001

Shower Chair Bench Heavy Duty



CH7955

Secure, Easy to Use – Our Toilet Safety Frame is approved the highest standard in medical Bathroom Toilet Rail, Tool-Free, Easy to assemble in 2 minutes

Barcode:6009701455618
Nappi Code:1128501001

Toilet Safety Rails with Anti Slip - Stand alone with backrest



CH7361LQH

This bath chair shower seat with arms assembles in just minutes.

- Total Height: 71 - 81cm
- Total length: 50cm
- Seat Width: 40cm
- Max weight: 130kg

Barcode:6009701455557
Nappi Code:1143639001

Shower Chair with Handles and No backrest

Bathroom Safety



Ch696

Material: Steel tube spraying
4 wheels with brakes, with footrest
Seat width: 46cm
Product size: 54*46*102-108cm
Net weight: 11kg

- Width between armrests 47cm
- Seat depth 46cm

Barcode :6009701451764
Nappi Code: 1103400001

Commode Mobile with footrest and castor lock





CH813

Serves as a safety frame – over an existing toilet
Easy to assemble
Comfortable seat
Weight capacity – 120kg

- Adjustable height from 56cm to 63cm floor to seat
- Adjustable height from 74cm to 84cm floor to armrest
- Width between armrests 47cm
- Seat depth 46cm

Barcode :6009701456042
Nappi Code:1143633001

Commode – With Drop Arm



CH799L

This high-quality shower chair is made of aluminum alloy.
The seat is contoured for comfort and stability and comes with drainage holes to prevent slippage.

- Seat Width: 68cm
- Seat Depth: 40cm
- Seat Height Min: 43cm
- Seat Height Max: 53cm
- Width: 80cm
- Depth: 50cm
- Weight: 5kg
- Maximum Weight Capacity: 130kg

Barcode:6009701451900
Nappi Code:1103404001

Shower Chair Bench With One Side Handle And Backrest



CH736LQ

Adjustable Height

- Clip on Armrests
- Max Weight: 130kg
- Seat Width: 40cm
- Total Height: 69-79cm
- Length Total:46cm
- Seat Height: 36 – 46cm

Barcode :6009701455533
Nappi Code:1103406001

Shower Chair with Back Support L



CH570

Both easy to fit and remove, it has adjustable feet that fit to the sides of the bath for added security.
Suitable for baths with an inside measurement of 40cm-63cm (16-25).

- Width 66cm
- Height 19cm
- Length 73cm

Barcode :6009701456059
Nappi Code:1143634001

Bath Board



Ch657

360 degree swivel to assist bathing the elderly, No fixtures or fitting required
Length 73cm

Barcode :6009701458800
Nappi Code:1166265001

Bath Chair - Swivel



Ch7945

Handrail for the bathroom can be positioned on any flat wall and comes in stainless steel, it is designed for elderly and disabled persons

The L-shaped / Shower Grab Rail provides a sturdy hand hold with anti slip rubber grip.
Net Weight: 5,82kg

Barcode :6009701456066
Nappi Code:1143635001

Hand Rail - Swivel

Bathroom Safety



CH5791 Double Stepper for Bath or Shower makes getting in and out of the bath tub easy and safe. Chrome plated steel pipe

Barcode :0710579687051
Nappi Code:

Double Stepper with Handle



CH736L

Easy to clean Height adjustable from 37-47 cm.
Screw in handles

- Unfolded width: 50cm
- Seat Width: 40cm
- Total Height: 71-81cm
- Total length: 47cm
- Max weight: 130kg
- Length Total: 46cm
- Seat Height: 36 – 46cm

Barcode :6009701451849
Nappi Code:1103407001

Shower Chair with Back Support S



Ch7956 Net Weight: 5.82kg

Grab Bar provides support for the elderly and helps people with limited mobility or physical disabilities as well

Barcode :0710579687068
Nappi Code:

Hand Rail L Shaped – Yellow



CHR02 Reacher

Barcode: 6009701459258

Grabber Tool



CH3641L

Swivel Shower Chair with handles and backrest
Pull lever to swivel and release to stop.
Lever on right hand side
Seat Width: 56cm
Seat Length: 52cm

Barcode:6009701451832
Nappi Code:1118678001

Swivel Shower Chair

Rollators



CH966LH

Lightweight
Foldable-convenient storage/transport
Padded Backrest- support
Bicycle style hand brakes-easy to use
Four wheels-aid mobility

- Unfolded Width: 62cm
- Seat Width: 45.5cm
- Total Height: 81-91cm
- Total Length: 62cm
- Seat Height: 56cm
- Wheel Dia: 6"
- Seat Depth: 32cm
- Backrest Height: 24cm
- Net Weight: 7.3kg

Barcode :6009701453430
Nappi Code:1103403001

6inch Rollator with 8 inch Solid Castors



CH965LH

Lightweight
Foldable-convenient storage/transport
Padded Backrest- support

- Unfolded Width: 62cm
- Seat Width: 45.5cm
- Total Height: 79-91cm
- Total Length: 69cm
- Seat Height: 56cm
- Wheel Dia: 8"
- Seat Depth: 32cm
- Backrest Height: 25cm
- Net Weight: 8kg

Barcode :6009701453423
Nappi Code:1103402001

6inch Rollator with 6 inch Solid Castors



CHWT04A

Foldable with high capacity shopping basket

- Frame Material: Steel tube
- Backrest: Sponge wrapped
- Cushion: Fabric Surface, sponge,wrapped
- Casters: 6pcs 6" universal wheel
- Load Capacity: 130kg
- Net Weight:11.4kg

Barcode :6009701455601
Nappi Code:1152070001

Walking Trolley



CH9181LH

Aluminum frame, Adjustable Handle height, Soft PVC Seat, Handle grips with Break system, 10"Front Castors

- Seat Width: 46cm
- Total Height: 81-93cm
- Total Length: 72cm
- Seat Height: 54cm
- Wheel Dia: 8"
- Weight Capacity: 130kg
- Net Weight: 8.2kg

Barcode :6009701453447
Nappi Code:1103408001

Rollator Supreme 4 wheel



CH914HW

Comfortable Seat: Compared to other walkers, this upright rollator features soft cushion, ideal seat for sitting. Dual handbrake system, ensure your safety when walking or sitting. The foot pedal, prevent from tipping over.

Barcode :6009701451092
Nappi Code:1120528001

4 Wheel Rollator with Armrest

Rollators



CHHWA61

Upgraded rehabilitation walking aid, the width of the double armpit can be adjusted, which is suitable for people of different body types. 30-42 cm, the length of the strap is adjustable, with a brake switch anti-rollover design;

Barcode :6009701452631

Hemiplegic Walker - Adult



CHHWC61

Upgraded rehabilitation walking aid, the width of the double armpit can be adjusted, which is suitable for people of different body types. 30-42 cm, the length of the strap is adjustable, with a brake switch anti-rollover design;

Barcode :6009701452587

Hemiplegic Walker - Child



CH969H

Product/Packaging Information
 Nett Weight: 7.4 kg
 Product Dimensions: 60cm(L) x 21cm(W) x 70cm(H)
 Shipping Weight: 8.0 kg
 Shipping Dimensions: 61cm(L) x 22cm(W) x 85cm(H)

Barcode :6009701451245
 Nappi Code:1106458001

3 Wheel Rollator with Bag and Basket

Aluminum Frame with adjustable handle and height. Soft PVC seat and handle grips, with break system. Detachable Backrest with Bag

- Seat Width: 46cm
- Total height: 54cm
- Length: 67cm
- Max Weight: 130kg
- Nett Weight: 7.3kg

Barcode:6009701453454
 Nappi Code:1103409001

Rollator 4 Wheel with bag and detachable back

CH9188LH



Walkers



CH961L

A Lightweight, Aluminium walker, provides a stable support for walking when more security is required than crutches. This version is equipped with a fabric seat.

- Unfolded Width: 55cm
- Seat Width: 46cm
- Total Height: 75-92cm
- Weight Capacity: 130kg
- Nett Weight: 3.4kg

Barcode :6009701451863
Nappi Code:1103465001

Walker with seat



CH9637L

A Lightweight, aluminum walker with 4 wheels provides a stable support for walking when more security is required than crutches.

- Unfolded Width: 55cm
- Total Height: 85-110cm
- Wheel Dia: 8"
- Load Capacity: 130kg
- Nett Weight: 4kg

Barcode :6009701451795
Nappi Code:1103461001

Walker with 4 wheels



CH913L

Folding Walking Frame Weighs less than 3 kg. Height adjustable

- Unfolded Width: 57cm
- Seat Width: 46cm
- Total Height: 56-73.5
- Nett Weight: 2.4kg
- Maximum Weight Capacity: 130kg

Barcode:6009701451627
Nappi Code:1103436001

Standard Walker



CHMK32

The Knee Scooter has large 20 cm wheels to roll easily over small bumps. It features a comfortable padded knee rest, no-tool height adjustment and brakes for safety. It also features a basket so you don't need to carry items in your hands.

Barcode:6009701458480
Nappi Code:1166309001

Knee Walker Foldable



CH912L

A Lightweight, aluminum walker, provides a stable support for walking when more security is required than crutches. The front wheels allow you to move the walker without lifting it.

- Unfolded Width: 63.5cm
- Total Height: 78-96cm
- Wheel Dia: 4inches
- Weight Capacity: 130kg
- Nett Weight: 2.75kg

Barcode :6009701455601
Nappi Code:1103975001

Walker with Wheels Foldable



CH50

Barcode : N/A
Nappi Code:6009701456189

Bottle Holder for Walker



CH9145HW

Aluminum frame and height-adjustable. Height adjustable Strong frame for outside moving Push handle provided with plastic grip

- Unfolded Width: 53cm
- Total Height: 48-58cm
- Net Weight: 2.8kg
- Weight Capacity: 100kg

Barcode :6009701456196
Nappi Code:1143640001

Walker - Non slip Hand Grip Paediatric

Wheelchairs



CH6926

A portable and discreet commode chair with a padded seat, backrest and armrests.

- Unfolded Width: 55cm
- Seat Width: 44cm
- Total Height: 95cm
- Total Length: 90cm
- Seat Height: 53cm
- Weight Capacity: 130kg
- Nett Weight: 12.5kg

Barcode: 6009701453683
Nappi Code: 1103445001

Wheelchair Commode with Foot Rest



CH809

Chromed Steel Frame, Fixed arm & Footrest, Excellent quality

- Seat Width: 46cm
- Total Height: 87cm
- Seat Height: 48cm
- Rear Wheel Dia: 24"
- Front Wheel Dia: 8"
- Backrest Height: 39cm
- Maximum Weight Capacity: 130kg
- Nett Weight: 16.6kg

Barcode: 6009701451238
Nappi Code: 1106454001

Standard Wheelchair




CHNM10B

Size 90*55*87cm
Wheel: Front 7inch
Rear: 22inch Aluminum pneumatic tire
Frame: Steel
Seat Width: 38cm
Seat Depth: 43cm
Paddles: Steel
Load Capacity: 120kg
Packaging: 90*26*87cm
Net Weight: 14kg
Gross Weight: 16kg

Barcode: 0710579687273
Nappi Code: N/A

Wheelchair Children- Manual Yellow



CH807

Powder Coated steel frame
Fixed Armrest
Folded Width: 21cm
Seat Width: 46cm
Total Height: 87cm
Seat Height: 48cm
Rear Wheel Dia: 24inch
Front Wheel Dia: 8inch
Total Length: 104cm
Weight Capacity: 130kg

Barcode : 6009701456028
Nappi Code: 1143641001

Wheelchair Aluminium Black Canvas



CH951B-51

Frame: Steel / Foldable, Upholstery: Nylon
Armrest: Removable / Desk
Footrest: Removable ABS Footplate
Caster:

- Seat Width: 51cm
- Total Height: 94cm
- Seat Height Max: 52cm
- Rear Wheel Dia: 24"
- Front Wheel Dia: 8"
- Total Length: 115cm
- Backrest Height: 40cm
- Maximum Weight Capacity: 170kg

Barcode: 6009701451665
Nappi Code: 1119974001

Wheelchair Deluxe Heavy Duty



CH976ABJ

Wheelchairs for adults and seniors make independent mobility enjoyable via specialized features

- Seat Width: 42.5cm
- Total Height: 101cm
- Seat Height: 52cm
- Rear Wheel Dia: 8"
- Front Wheel Dia: 8"
- Weight: 13kg
- Maximum Weight Capacity: 130kg

Barcode: 6009701451054
Nappi Code: 1120526001

Transit Wheelchair - Red

Wheelchairs

CHDX01



A portable and discreet commode chair with a padded seat, backrest and armrests.

- Motor Power 24V, 250w
- Controller: 24V/45A
- Front Wheel: 9 inch solid tyre
- Load Capacity: 150kg
- Climbing height: 12 inch
- Endurance: 18-20km
- Speed: 0-8km/h
- Shock Damping: yes rear damper
- Light: Front LED
- Seat: PU

Barcode:6009701456035
Nappi Code:1143642001

Scooter Electric - with 360 degree moving seat

CHDX02



Flip Back, with adjustable armrests, swivel seat & horn.

- Motor Power: 24V, 250w
- Controller: 24V/45A
- Battery:24V 12AH*2pcs lead acid battery
- Front Wheel: 9 solid tyre
- Product Size: 112*52*90cm
- Load Capacity: 150kg
- Climbing height: 12"
- Endurance: 25-28km
- Net weight:55kg(including battery)
- Shock Damping: yes rear damper

Barcode: 6009701457315
Nappi Code:

Scooter Electric - with detachable rear

CH505



This is suitable for wheelchairs with a distance of 19.5 to 20.5 inches between the two armrests. But remember, for portable wheelchairs, manual wheelchairs; not for electric wheelchairs, power wheelchairs, round armless wheelchairs.

Material: ABS plastic

Barcode:6009701456165
Nappi Code:N/A

Wheelchair Tray

CHSP02



9 inch Solid Tyre replacement for Scooter

Barcode:6009701459296

Tyre - Spare for Scooter

CHM04



Powder Coated Frame,Fixed Arm & Footrest, Foldable, Solid Rear Wheel

- Frame Material: Carbon Steel
- Unfolded Size: 90*60*89cm
- Folded Size: 74*27*69cm
- Seat Width: 45cm
- Front wheel: 6
- Rear Wheel: 20 Solid Tyre
- Safety Load: 110kg
- Handrail: Soft PU
- Cushion Material: Oxford Fabric
- Cushion Colour: Red

Barcode:6009701456394
Nappi Code:N/A

Wheelchair Children's Standard Red Foldable

CHBT01



Non-Cadmium Design, Environment-friendly: Chilwee Battery has adopted internationally leading technology

Barcode:6009701456165

Battery for Scooter - Lead Acid 24V

CHM17



High Backrest reclining by 90degrees to 180degrees
Detachable flip up PVC, armrests combined with detachable elevating footrests

- Seat Width: 45cm
- Total Height: 128cm
- Seat Height Max: 50cm
- Rear Wheel Dia: 24"
- Front Wheel Dia: 7"
- Total Length: 119cm
- Maximum Weight Capacity: 130kg

Barcode:6009701456387
Nappi Code:

Wheelchair with high back Full recline with Commode

CHM13



Drop Down Handles
Steel frame Drop Back Handles for Compact Storage.

- Seat Width:45cm
- Front wheel: 7" Solid Tyre
- Rear Wheel: 24" Solid Tyre
- Seat Height: 45cm
- Max Weight: 130kg
- Handrail: Soft PU
- Break Method: Double Break (front and rear)

Barcode:6009701456318
Nappi Code:1143643001

Wheelchair with drop down handles

Wheelchairs

CH02B



Unfold Size: 96*64*89cm
 Folded Size: 96*27*89cm
 Seat Width: 45cm
 Front wheel: 6" Solid tyre
 Rear wheel: 20" solid tyre spoked wheel
 Safety Load: 135kg
 Backrest: collapsible
 Handrail: ABS
 Foot Rest: Aluminium alloy, height adjustable
 Break method: four wheel brakes
 Cushion material: mesh cloth
 Cushion colour: Orange or Blue
 Net Weight: 13kg

Barcode: 6009701459272
 Nappi Code:

Wheelchair Aluminium Blue/Orange and Black - 45cm

CHM01



Frame Material: Carbon Steel
 Unfolded Size: 82*61*83cm
 Folded Size: 74*28*69cm
 Seat Width: 40cm
 Front wheel: 6
 Rear Wheel: 16inch Solid Tyre
 Safety Load: 110kg
 Handrail: Plastic

Barcode: 6009701458022
 Nappi Code: N/A

Wheelchair Childrens Standard
 Black Foldable 41cm

CHE05



- Frame Material: Carbon steel
- Motors: 250W*2 Electronic Brake System
- Charging voltage: 110-220V
- Charging Time: 6-7 hours
- Edurance: 15-20km
- Max Speed: 6km/h
- Front Wheel: 10" Solid Tyre
- Rear wheel: 16" pneumatic tyre aluminium
- Seat Width: 45cm
- Safety Load: 120 - 140kg
- Product Size: 120x62x128cm
- Battery: 24V, 12AH, lead acid battery
- Net Weight: 47kg (including battery)
- Shock absorption: Front and rear

Barcode: 6009701456165
 Nappi Code: N/A

Luxury Electric Wheelchair

CH9602



Versatile exercise machine
 Helps to improve the blood circulation
 Adjustable level of resistance
 Stylish design
 Technical Specifications:

Barcode: 6009701456141
 Nappi Code:

Pedal Exerciser

Wheelchairs

CHB16



- Frame Material: aluminum alloy frame
- Motors: 250Wx2 Electromagnetic braking system
- Charging voltage: 110-220V
- Endurance: 10-15km (Depending on road condition and weight)
- Max Speed: 6km/h
- Climbing ability: <13degrees
- Front wheel: 12inch Pneumatic tyre
- Seat Width: 47cm
- Safety Load: 100kgs
- Product size: 103x65x90cm

Battery: 24v, 12AH Lead acid battery

Net Weight: 32kg (including battery)

Shock absorbent front damper

Barcode: 0710579687143
Nappi Code:N/A

Electric Wheelchair – Folding Lightweight

CHWR5



Wheelchair Ramp - 4 Feet 1.2 Meters

Barcode: 0710579687433
Nappi Code:N/A

Wheelchair Ramp - 5 Feet 1.5 Meters

CHWR4



Wheelchair Ramp - 4 Feet 1.2 Meters

Barcode: 0710579687440
Nappi Code:N/A

Wheelchair Ramp - 4 Feet 1.2 Meters

CHE19



- Frame Material: Aluminum alloy frame
- Motors: 250W*2 Electromagnetic braking system
- Charging Voltage: 110-220V
- Endurance: 10-15km (Depending on road condition and weight)
- Max Speed: 6km/h
- Climbing ability: < 13
- Front Wheel: 12 pneumatic tyre
- Seat width: 47cm
- Safety Load : 100 – 140kg
- Product size: 104x60x102cm

Battery: 24V, 12AH Lead Acid Battery

Recline available: Yes

Barcode: 0710579687150
Nappi Code:N/A

Electric Wheelchair with Headrest

CH976ABB



Wheelchairs for adults and seniors make independent mobility enjoyable via specialized features

- Seat Width: 42.5cm
- Total Height: 101cm
- Seat Height : 52cm
- Real Wheel Dia: 8"
- Front Wheel Dia: 8"
- Weight: 13kg
- Maximum Weight Capacity: 130kg

Barcode:6009701459616
Nappi Code:

Transit Wheelchair - Black

CH507



The non slip surface of the wheelchair tray for adults ensures that items placed on it stay in place, preventing accidents and spills.

Barcode: 0710579687709

Wheelchair Tray – Black

CH506



Durable Plastic Tray: Crafted from high-quality White ABS plastic The tray table offers durability without adding unnecessary weight to your wheelchair.

Barcode:0710579687716

Wheelchairs

CH01B



Aluminum alloy meta frame
 Product Size: 84x48x89cm
 Seat Width: 40cm
 Front Wheels: 6 inch solid PVC Material
 Rear Wheels: 8 inch solid PVC Material
 Safe Loading Capacity 100kg
 NW: 7.5kg

Barcode: 0710579687143
 Nappi Code:N/A

Wheelchair Aluminium Transit

CHDX03



Anti Tip Wheels
 8-14 hours daily use
Battery: Two Deep Cycle batteries / AGM and Gel Cell
Maximum Slope: 10 degrees
Max Climbing: 12 Degrees
Maximum Obstacle Climbing Ability: 102cm
Overall Length: 102cm
Overall Width: 51cm
Speed: Variable up to 8.00km/h
Weight Capacity: 120kg
Seat dims: 44cm width
 1 step folding button
Net Weight: 34kg (including battery)
 Shock absorbing front damper

Barcode:
 Nappi Code:N/A

Folding Scooter Z ONE 4 Wheel

CHDX04



- Front lighting
- Adjustable Shock-Absorbing Tiller
- Large visible display screen
- 10Ah Airline friendly lithium-Ion battery for airtravel
- Front wheel with 8 PU foam solid tire Rear wheel with 9 PU foam solid tire
- Optional storage basket
- Removable and flip-up armrests
- Ergonomically designed seating position
- Highly visible rear lighting
- Type-C and USB port for charging mobile devices
- Arm rests
- Safety belt
- Under-Seat Basket

Barcode:
 Nappi Code:N/A

Scooter Compact Folding - RooBooter Q40

Beds and Mattresses




CH030-125

Patient Examination Bed, is used by Doctors, Care Givers, First Aiders, in the examination and small treatments of patients. Foldable/Black Comes with a easy to carry bag

- Net Weight: 2.35kg
- Gross Weight: 2.44kg

Barcode:6009701453577
Nappi Code: N/A

Examination / **Massage Bed - Black**



CHMASP1

Stainless steel base with a purple cushion

Dimensions: 186cm (L) x 68cm (W) x 58-85cm (H)

Barcode:6009701458817
Nappi Code:N/A

Massage Bed - With backlift function Purple





CHMASF2

Stainless steel base with a black cushion

Dimensions: 186cm (L) x 68cm (W) x 58-85cm (H)

Barcode:6009701453577
Nappi Code: N/A

Massage Bed - Black without backlift function



CH2561

- Bedboard : overall punching, carbon steel
- Head&foot board: ABS with fender wheels
- Siderails: 6 bars aluminum alloy with metal lock
- Casters: central locking system double sides casters
- Crank: ABS foldable
- Mattress: 6cm thickness with waterproof cover
- Infusion pole: stainless steel with 4 hooks
- **Two Functions:**
- Back adjustment 0°~80° 2. leg adjustment 0°~45°
- **Outer Size:**
- 2060*1100*500mm
- Load Capacity: 250kg
- ABS with fender wheels, 1. 6 bars, aluminum alloy with metal lock
- 40cm height central locking system, double sides, casters, ABS, foldable.
- IV: stainless steel with 4 hooks

Barcode:6009701456080
Nappi Code:N/A

Two-crank Abs Hospital Bed With 4 Section Mattress



CH633

- Bedboard: overall punching, carbon steel
- Head&foot board: PP/ABS
- Siderails: 6 bars aluminium alloy with plastic lock,40cm height
- Casters: 4pcs - ABS noiseless caster with single lock
- Crank: 3pcs - ABS foldable
- Mattress: 6cm thickness with waterproof cover, palm and sponge inner
- Infusion pole: 1pc,stainless steel with 4 hooks, height adjustable
- **Three Functions:**
- back adjustment 0°~80°
- Leg adjustment 0°~45°
- Height adjustment 400°~700mm
- Outer Size: 2060x1000x400x700mm
- Load Capacity: 250kg
- Bed Packing size:203x 96

Barcode:6009701453829
Nappi Code:N/A

Three Crank Abs Hospital Bed With 4 Section Mattress

Beds and Mattresses



CHS103

This is a attendant bed for a hospital cabin. Fully MS steel frame and high quality mattress. Foldable bed come chair for anywhere use. Bed turns into a chair. It has a powder coating paint and wheel for easy movement.

Barcode:6009701456097
Nappi Code:N/A

Sleeping Chair



CHM09A

- Bedboard: banded, carbon steel
- Head and foot board: PP
- Siderails: 5 bars Aluminum Alloy with plastic lock
- Casters: ABS noiseless caster with single lock
- Crank: electroplate metal foldable crank
- Mattress: 6cm thickness with waterproof cover
- Infusion pole: stainless steel with 4 hooks
- Dining table: ABS
- Packing size: 203*95*19cm/1pc

Barcode:6009701456097
Nappi Code:N/A

**One Crank ABS Hospital Bed
Mattress - Manual**



CHCSST07

Stainless Steel Frame
4cm thick mattress

- Single Crank
- Back lifted
- Size: 1900x600x700mm

Packaging SizeL 192x62x20cm

Barcode:6009701451825
Nappi Code:N/A

Examination Bed 1 Crank



CA-H05B

Collapsible Retractable Aluminum
Side Rails /Hospital bed Furniture

Product

Material Aluminum alloy
Specification 1200L X 350H mm
One-touch type handle operation
Aluminum alloy folding guardrail

Barcode:6009701453577
Nappi Code: N/A

Side Rails - Crank Bed



CH2562

Bedboard: overall punching carbon steel
Head and Foot board: ABS and stainless Steel composite detachable with fender wheels
Siderails: 5 Bar stainless steel and metal lock 40cm height
Castors: 4pcs ABS High quality noiseless caster with single lock
Crank: 2Pcs electroplate metal foldable crank
Infusion Pole: 1pc stainless steel with 4 hooks

- **Two Functions:**
- Back adjustment 0°-80° 2. leg adjustment 0°-45°
- **Outer Size:**
- 2060*1100*500mm
- Load Capacity: 250kg
- ABS with fender wheels, 1. 6 bars, aluminum alloy with metal lock
- 40cm height central locking system, double sides, casters, ABS, foldable.
- Mattress 6cm thickness with waterproof cover
- IV: stainless steel with 4 hooks

Barcode:6009701458503

Nappi Code:N/A

**Two-crank Abs and Steel Hospital Bed With
4 Section Eggbox Mattress - Manual**



GYNBEDO1CH

- Carbon Steel Frame
- 1.2mm diameter thickness
- Function: Back and leg adjustable

AVAILABLE IN PINK

Barcode:6009701456301 Pink Barcode:0710579687754
Nappi Code: N/A

Gynaecological Examination Bed

Beds and Mattresses



CHBD3

High strength carbon steel pipe + PU armrest cover + PU mesh storage bag
 Armrest Length: 73 - 93cm
 Height: 75 - 93cm
 Height Adjustable in 15 steps
 Size 37,5cm
 The 2 sides of the armrest are adjustable in 2 steps and the adjustable range on one side is 10cm.
 Gross Weight: 5kg

Barcode: 6009701459241
 Nappi Code: N/A

Bedside Handrail



CHEXCB

- Black cover
- Stainless steel frame
- Product size: 180*60*75cm
- Net weight: 42kg

Barcode: 6009701459128
 Nappi Code: N/A

Examination Couch Standard - Black



CHE03

With 8cm eggbox sponge+PU cover
 Electric 3 functions Parameter

1. Bedboard: overall punching, carbon steel material
2. Motor: 3pcs, DC linear motor with remote controller
3. Head&foot board: ABS with fender wheels
4. Siderails: 6 bars aluminium alloy with plastic lock, 40cm height
5. Casters: 5 inch diameter, "DABEN" brand ABS high quality noiseless casters
6. Infusion pole: 1pc stainless steel with 4 hooks
7. Urine bag hooks

Mattress: 190*90*40cm(5pcs/bag)
 Total Bags: 5bag
 Total Volume: 3.3cbm
 Total Weight: 145 Kg

Barcode: 6009701458503
 Nappi Code: N/A

Three Function Electric Bed with 8 cm Eggbox Mattress - Manual



CHPL1

Material: Carbon Steel
 Product Size: 1050(L)x1000(W)mm
 Height Adjustable: 01 - 1850mm

Barcode: 6009701456301
 Nappi Code: N/A

Patient Lifter - Battery Operated

Beds and Mattresses

CHN23A Standard Hospital Bed



Barcode: 6009701451023
Nappi Code: N/A

Hospital Bed Standard

PTLFCH



• Freestanding lifting pole for patient 0808195 to raise the user from supine to sitting position, the base is placed on the ground and the hanger is attached to it.

Dimensions of 79X170X75CM
Suitable for a maximum weight of 110kg.

Barcode: 6009701459814
Nappi Code: N/A

Patient Lifting Pole - Free Standing

CHPETWS




Product Size: 2000 x 600 x 650mm
Material: Medical antibacterial board
High Quality, soft material leather, skin friendly comfort, breathable and non-allergic

Barcode: 6009701459791
Nappi Code: N/A

Pediatric Exam Table with stairs

BACAR01CH



Bassinet Trolley Complete
Material: Stainless Steel
Includes Iv Pole

AVAILABLE IN BLUE

Barcode: 6009701459326
Nappi Code: N/A

Stainless Steel Baby Carriage

CHS211 Bedboard: banded carbon steel
Product size: 160x80x63cm



Barcode: 6009701452471
Nappi Code: N/A

Child Cot Bed with Crank

Beds and Mattresses



CH210

Child Cot Bed Mobile with drop down cotsides
Child Cot Bed with drop down cotsides Hospital Bed for children in a Paediatric Ward.
Safety Guard Rails surround the bed, and can be dropped down completely. This allows the safety of the child patient.

- Material: Carbon Steel
- Drop Down cot sides on both sides
- Dimensions: L 80x50x55cm

Barcode:6009701450934
Nappi Code:1118294001

Child Cot Bed with Mattress



CHEGBX

Egg Box Mattress 4 Part 10cm
for 2 and 3 crank hospital beds
10cm thickness
PU Cover
Net weight: 5kg

Barcode:6009701456936
Nappi Code: N/A

Egg Box Mattress 4 Part 10cm



CH10WP

TPU Mattress (4 sections)
for 2 and 3 crank hospital beds
10cm thickness 33A High density sponge
TPU Material waterproof cover
Net weight: 6kg

Barcode: N/A
Nappi Code: N/A

Mattress 10cm 4 Section Waterproof



CH592

Bubble mattress is designed for the treatment and prevention of pressure sore stage 1 in short-term and domiciliary therapies. The mattress is PVC is composed of 130 7cm bubble cells, particularly comfortable, it is fixed to the bed by extra flaps on both top and bottom side. Very easy to install, maintain and use thanks to the possibility of adjust comfort range.

Odourless Medical-Grade PVC
Japan One-Time Forming HF Device
High Quality Rubber Tube
Adjustable hangers to suit any kind of bed frame
Max User weight: 135kg
Dimensions without flaps (L x W x H): 200X 90 X 7cm
Flap length: 50 + 50cm
Material: Odourless Medical-Grade PVC.
Cycle Time: 6 Min

Barcode:6009701450934
Nappi Code:1118294001

Bubble Mattress with Pump



CHS303

The medical anti-bedsore mattress consists of a number of bars or spherical chambers.
It is arranged into a single double exhaust chamber in a double-row arrangement, which can change the pressure point of the body continuously, and keep the blood circulation of the patients in normal circulation.
There are regular micropores in the air chamber, which can be pumped out of the micropores to keep the skin dry, thus preventing bedsores.
Air Mattress size: 1950 x 900

Main Technical Indicators
Charge and release cycle: 5min for each aeration, 5 min for air release, 10min for one cycle, plus or minus 5% of the gas cycle
Maximum inflatable pressure: the air charging pump can be adjusted by inflating pressure with maximum pressure greater than 12kpa:
Maximum bearing weight: 135kg
• Gas pump: no more than 20VA

Barcode:6009701455649
Nappi Code:1133707001

Anti Bedsore Electric Mattress

Beds and Mattresses

CHBLFT



Bedboard size* 1900mm x 900mm
 2cmx3cmx1.5mm pipe
 Mattress: 7cm thickness foam
 Functions: Backboard lift, leg board lift

Barcode:0710579687730

Electric Back Lifter

CH2560



- Bedboard: Banded carbon steel
- Head and Foot board: ABS detachable
- Siderails: 5 Bar aluminium alloy with plastic lock 40cm height
- Castors: 4pcs ABS High quality noiseless caster with single lock
- Crank: 2Pcs ABS Foldable
- Infusion Pole: 1pc stainless steel with 4 hooks
- Urine bag hooks

Barcode:6009701458497

Nappi Code: N/A

Two-crank Abs Hospital Bed With 4 Section Eggbox Mattress

CHLG10



- Black cover
- Stainless steel frame
- Product size: 180*60*75cm
- Net weight: 42kg

Barcode: 0710579687747

Nappi Code: N/A

Leg Pillow

CHHA11



Heavy Duty
 Crank lever for adjusting
 Dimensions: 190cm x 82cm x 54cm –
 minimum height
 190cm x 82x 76cm maximum height
 Features: Face hole, crank
 Colour: Blue

Barcode:0710579687723

Nappi Code: N/A

Examination Bed Heavy Duty Height Adjustable

CHP300



Maximum inflatable pressure: the air charging pump can be adjusted by inflating pressure with maximum pressure greater than 12kpa: Air mattress under maximum working load, the maximum pressure is greater the 4kpa
 Maximum bearing weight: 135kg
 Gas pump: no more than 20VA

Barcode: 6009701459784

Nappi Code: N/A

Anti Bedsore Pillow for wheelchair

Crank Bed
 Head Board and Foot Board



Crank Bed Foot Board
 ABSCRFB01

Crank Bed Head Board
 CRBHB01

Crank Bed - Foot and Head Board

Stretchers

CH00202



Product size:195*50*18cm
 Fold size: 100*12*35cm
 Frame material: high strength stainless steel material and Oxford cloth material
 This stretcher has legs and can be folded.
 With light weight, small size,easy to carry, reliable quality, safe use, easy to sterilize.
 Color options: blue
 Packing size:105*40*27cm(5pcs)
 Net weight: 5.5kg

Barcode:6009701453577
 Nappi Code: N/A

Stretcher - Blue 2 Fold with wheels

CH00303



4 Fold aluminum stretcher
 Includes carrier bag
 Orange

Barcode:6009701455632
 Nappi Code:1128499001

Stretcher - Orange 4 Fold with bag

TRAN03CH



This board is used to help a person easily and safely transfer from a wheelchair, bed, chair, sofa, commode or other sitting position while minimizing user and caregiver strain or injury

Barcode:6009701458510
 Nappi Code:N/A

Transfer Board

CH035




The sturdy aluminum frame is strong enough to hold up to 160kgs
 Three adjustable straps with quick-release buckles ensure the patient's comfort and safety..
 Non-slip grips on the rear and front handles help rescuers maneuver the Aluminum Stair Chair easily

- Product size: 90*51*91cm
- Fold Size: 90*17*59cm
- Load Capacity: 160kg
- Frame Material: high-strength aluminum alloy
- Adopts folding structure
- Stretcher equipped with 2 or 4 wheels
- Two or four foldable handles are provided at the back end.

Barcode:6009701457940
 Nappi Code:N/A

Stretcher
 Stair Chair with 2 wheels

CHPM8



Used for the transportation of casualties inflicted with spinal injuries.

- Product Size: 183*43*5cm
- Load Capacity: 160kg
- Frame Material: high-strength engineering plastics
- Spine boards are commonly used to safely transport severely injured and immobile people.
- If floats on water, allows X-Rays and can also be used with a head retainer.
- This spine board is durable, age-resistant, lightweight and easy to store.
- It is essential for emergency equipment in medical centers.

Barcode:6009701457957
 Nappi Code:N/A

Spine Board

CHOP02



Used specifically for moving injured people. It is most frequently used to lift people who may have a spinal cord injury from the ground, either due to unconsciousness or in order to maintain stability in the case of trauma.

- Product Size: 208*44*7cm
- Fold Size: 120*44*9cm
- Load Capacity:160kg
- Frame Material: high strength aluminium alloy
- Both ends are equipped with clutch device which can be separated into left and right side switch and easier to operate.

Barcode:6009701457964
 Nappi Code:N/A

Scoop Stretcher
 Foldable extension type

Stretchers

CH00101



Double folded aluminum stretcher

Barcode:6009701458459
Nappi Code:N/A

Stretcher
Orange 2 Fold with Bag

CHS10



A detachable stretcher for ambulances, the bed surface can be quickly separated, suitable for transportation inside the hospital and emergency center transportation and outdoor rescue. With the function of height lifting and back adjustment, it can well meet the needs of various scenes.

Main features

- This multi functional automatic loading stretcher has varied heights, it's made of high-strength aluminum alloy thick-wall pipes.
- Support for IV stand included
- Aluminum alloy frame, easy cleaning, and sterilization.

Barcode:6009701458848
Nappi Code:N/A

Stretcher - Emergency/Ambulance

CHPT22



The product was made from high quality stainless steel;
Beautiful stainless steel foldable guardrails;
Four general casters;
Good quality
Size:1900x550x700mm

Barcode:6009701451016
Nappi Code:N/A

Stainless Steel Flat Emergency Patient Stretcher with Backrest Adjustable

CH036



The sturdy aluminum frame is strong enough to hold up to 160kgs. Three adjustable straps with quick-release buckles ensure the patient's comfort and safety. Non-slip grips on the rear and front handles help rescuers maneuver the Aluminum Stair Chair easily.

- Product size: 90*51*91cm
- Fold Size: 90*17*59cm
- Load Capacity: 160kg
- Frame Material: high-strength aluminum alloy
- Adopts folding structure
- Stretcher equipped with 2 or 4 wheels
- Two or four foldable handles are provided at the back end.

Barcode:6009701458831
Nappi Code:N/A

Stretcher - Stair Chair with 4 Wheels Black

EMRSCCH



First Aid Rescue Sheet
Protection against cold, heat and humidity
Individually packed in a mini grip bag
Packing size:105*40*27cm(5pcs)
Net weight: 5.5kg

Barcode:0710579687556
Nappi Code: N/A

Emergency Rescue Blanket - Singles

Stretchers

CHS11



Soft Stretcher - Orange with a handbag
 Product Size: 212x75cm
 Net Weight: 9kg
 Gross Weight: 10kg

Barcode:6009701459647
Nappi Code: N/A

Soft Stretcher

CHSK21



Splint Kit / Set
 EVA Material
 6pcs/set
 Gross weight: 2,3kg

Barcode:6009701459630
Nappi Code:

Splint Kit / Set

CHCT103



It is a manual transfer stretcher.
 Sleep surface : 1900x610mm
 HI / LO: 875-570mm,
 Backrest Lifts angle: 0-60
 Central control lock system castors
 Net weight: 84kg

Barcode:6009701459333
Nappi Code:N/A

Transfer Stretcher

CHPTC1




Foot control, by hydraulic adjustment
 Underarm support design to lift the user safely.
 Material: carbon steel spraying
 Product size: 63x55x118-138cm
 Seat height adjustment: 42x50cm
 Toilet hole size: 21x21cm
 3 inch noiseless wheel

Barcode: 6009701459623
Nappi Code:N/A

Patient Transfer Chair

CHPM9



Our spine boards are designed with safe belts made of HDPE that is pollution-free and remarkably wear resistant.


The patient handling device is floatable, so it can be used for water rescue.

Because of its translucency, it allows X-rays to pass through.

Barcode:6009701459852
Nappi Code:N/A

Spine Board with Spider Strap and head holder + Bag

CHPTST1



Barcode: 0710579687327
Nappi Code:N/A

Standing Transfer Chair

CHSK20



Complete splint system that contains arm splint, adult & child; leg splint, adult & child; arm/shoulder immobilizer

In-line fiberglass support for transport and stabilization

Barcode:
Nappi Code:N/A

Splint Kit - TPU Orange

Medical Furniture

CHB075-34

- Medical screen 4 fold
- 4 fold, with 7 wheels
- Stainless Steel frame
- Waterproof woven fabric
- Product size: 2000*1800mm
- Gross weight: 20kg



Barcode:6009701453584
Nappi Code:N/A

Bed Screen 4 Fold

CHB075-33

- 3 fold, with 6 wheels
- Stainless Steel frame
- Waterproof woven fabric
- Product size: 1500*1800mm*
- Gross weight: 15kg



Barcode: 6009701456264
Nappi Code: N/A

Bed Screen 3 Fold

CHZ720

Thick Pad- Brown egg box finish
It help you sitting comfortable even after many hours.
Rolling Smoothly
5 Wheels
Easy to assemble/disassemble.
Height Adjustable




Barcode:6009701458534
Nappi Code:

Salon Stool

CHPSST25

- Material: stainless steel frame
- PU seat cover +sponge inner
- Seat height adjustment: 47-62cm
- Seat diameter: 34cm
- Back height: 20cm
- Net weight: 3.6kg
- Material: ABS plastic



Barcode:6009701458510
Nappi Code:N/A

Chair Nurses Metal Black with wheels

CHZF25

- Thick Pad—Seating diameter:15.75 inches,Thickness:4.7 inches;
- Rolling Smoothly—360° Swivel caster and ultra-quiet wheels allow for easy movement.
- Soft polyurethane wheels can prevent scratches on your floor.
- Easy to assemble/disassemble.
- Height Adjustable



Barcode:6009701456134
Nappi Code:N/A

Doctors Chair

EXLM23CH Examination Lamp



Barcode:6009701458527
Nappi Code:

Examination Lamp - LED

CHZ101

- Infusion Stand with handle
- Stainless steel pipe
- Height adjustable
- 4 wheels with brake
- 4 hooks



Barcode:6009701451764
Nappi Code:

IV Stand with Handle

CHG01

Frame Material: Steel Tube spraying
92x65x34cm

Adjustable Backrest
Size (L) - 1000 X (W) - 660 x (H) - 14100 mm



Barcode:6009701457339
Nappi Code:N/A

IV/Infusion Reclining Chair

Medical Furniture




CHRD1

Full Body Manikin set
 1 x full body manikin
 1 x Monitor
 1 x Portable plastic backing box
 1 x User Manual
 1 x CPR operation Pad
 50 x CPR face shield sheet
 4 x Exchangeable lung sets
 2 x Rolls of Thermal printer paper
 1 x Data Cable
 1 x Adapter

Barcode: 6009701459340
Nappi Code: N/A

Full Body CPR Doll



CH00101

Packaging Dimensions: 800*800*200

5 wheel base
 Accessory basket
 Weight: 6kg
 Mounting Plate

Barcode: 6009701458619
Nappi Code: N/A

Patient Monitor Trolley - Hydraulic



CHZ721

- Thick Pad- black egg box finish
- It help you sitting comfortable even after many hours.
- Rolling Smoothly
- 5 Wheels
- Easy to assemble/disassemble.
- Height Adjustable

Barcode: 6009701459838
Nappi Code: N/A

Salon Stool - Black




CHZF26

Doctors Chair hydraulic control
 Frame material: Aluminium
 Seat height adjustment: 43-57cm
 Base Diameter: 5.6cm
 Gross weight: 7.5kg

Barcode: 6009701459821
Nappi Code: N/A

Doctors Chair with Backrest



CHEXBL2

Examination Lamp designed to revolutionize medical lighting with its advanced features and exceptional performance.

This lamp is meticulously engineered to enhance the examination experience for medical professionals while ensuring utmost patient comfort and safety.

Barcode:
Nappi Code: N/A

Examination Lamp – Black



CHEXWH2

Examination Lamp designed to revolutionize medical lighting with its advanced features and exceptional performance.

This lamp is meticulously engineered to enhance the examination experience for medical professionals while ensuring utmost patient comfort and safety.

Barcode:
Nappi Code: N/A

Examination Lamp - White

Medical Furniture

CHRDL



Full Body Manikin set
 1 x full body manikin
 1 x Monitor
 1 x Portable plastic backing box
 1 x User Manual
 1 x CPR operation Pad
 50 x CPR face shield sheet
 4 x Exchangeable lung sets
 2 x Rolls of Thermal printer paper
 1 x Data Cable
 1 x Adapter

Barcode: 0710579687778
 Nappi Code: N/A

Full Body CPR Doll - Child

INTRSSCH



As per the changing and ever rising needs of our customers, we are providing a comprehensive collection of SS Dressing Tray in multiple specifications. Used by nurses during operations and surgeries, these offered trays are light in weight and are easy to carry. Besides this, these are developed using sturdy and durable material under the direction of skilled professionals.

Barcode: 0710579687761
 Nappi Code: N/A

Instrument Tray Stainless Steel with Cover 24x16x5cm

CHMNLMP



Barcode: 0710579687334
 Nappi Code: N/A

Magnifying Lamp

CHMNLMP



Barcode: 0710579687303
 Nappi Code: N/A

Waiting Chair

Trolleys

CHC074



- Thick stainless steel, brushed board surface, not easy to scratch the hand, not easily deformed, anti-corrosion rust.
- Mute slide single drawer design, stainless steel, no smell, easy to clean, adjustable plastic dirt bucket design.
- Comfort handle, stainless steel guardrail, seamless welding, wear-resistant, impact resistant.

Barcode:6009701453737
Nappi Code: N/A

Medical Trolley Large Stainless Steel

CH5069S



- Stainless steel tubular frame
- 2 tier stainless steel trays
- Dismantled for easy transportation –
- Easy to assemble
- Suitable for carrying multipots with plinth and drip trays

Barcode:6009701452600
Nappi Code:N/A

Stainless Steel Trolley - 2 Layer

CHC071



Stainless Steel Trolley with bucket

86x44x72cm

Barcode:6009701453829
Nappi Code:N/A

Clinic Trolley Large

CHPS02A



Aluminum Alloy Frame, Wooden Board, Hydraulic Lift, 4 Piece Noiseless Wheels with brake
Product Size: 82*45*(84-107cm)
Packing size: 86x54x27cm

Barcode:6009701456103
Nappi Code:N/A

Overbed Table - ABS with Wooden Top

CH5662



Overbed Table
Manual Hi/Low Adjustment
Epoxy
Dimensions of Table:
Width - 38cm
Length - 81.5cm
Depth: 1.8cm

Barcode:6009701458442
Nappi Code:N/A

Overbed Table - Wooden Top Swivel

The desktop is made of ABS material

The air-controlled lifting device can adjust the height at will according to personal needs and preferences.

4 x 2.5inch casters with brakes, which can move the position at will and make the docking more stable.

Barcode:6009701458541
Nappi Code:N/A

Overbed Table - ABS

ABOT02CH



Trolleys



CHMT12

Product Size: 710x620x1010mm
Wooden Carton
ABS

Barcode:6009701457353
Nappi Code: N/A

Treatment Trolley with 3 tier worktables and wheels - Blue



CHMT13

Product size: 710*620*1010mm
Wooden carton*

Barcode:6009701457360
Nappi Code:N/A

Treatment Trolley with 3 tier worktables and wheels - Pink



CHMT10

Three layers
ABS Net Weight: 42kg
Gross Weight: 46kg
Packing Size 780*560*940mm

Barcode:6009701457346
Nappi Code:N/A

Treatment Trolley ABS



CHS402

Durable strong plastic cabinets used to store multiple items or clothes of a patient at a hospital or clinic.

Overall Size: 480*480*760mm
Material: ABS Engineering plastic
With towers on both sides, one dinning board, one drawer and one door.
Packing Size: 500*500*780mm
N.W.: 15kg
G.W.: 19kg

Barcode:6009701451771
Nappi Code:N/A

Bedside Locker ABS with wheels



CHS403

- ABS rescue vehicle/emergency vehicle parameters
- Suitable for medical staff to give first aid to patients
- Has double bottom injection molding process to form the table, recessed design can prevent items from falling. 304 stainless steel three sided fence, the table is equipped with transparent soft glass, a sharp tool box.
- Back: Cardiopulmonary resuscitation board, oxygen cylinder holder, power cord socket
- Left Side: Assistant working desk, one data box
- Right Side: two dirt buckets, four luxury silent casters, (two with brake function, two without brake) anti curling hair
- Car body specifications: 750*480*920mm

Barcode:6009701456110
Nappi Code:N/A

Emergency Trolley - 5 Drawer



CH276

Stainless steel 3 drawers trolley
Product Size: 66*44*86cm
Packing Size: 70*45*80cm
Net Weight: 20kg

Barcode:6009701457933
Nappi Code:N/A

Stainless Steel Trolley - 3 Layer

Trolleys



CHLT1

Material: Stainless steel
 Color options: blue/yellow
 Product size: 900x500x760mm
 4 wheels with brakes
 Net Weight: 13kg

Barcode:6009701459104
 Nappi Code:N/A

Laundry Trolley



CHSSW10

Detachable type
 Material: stainless steel
 2 shelves, 1 big drawer, 4 wheels with brakes
 Product size: 66*44*86cm
 Net weight: 12kg

Barcode:
 Nappi Code:N/A

Trolley Stainless Steel with wheels - Wooden Top



ANETRCH1

ABS Trolley
 Product Size

Barcode:0710579687334
 Nappi Code:N/A

Anesthesia Trolley



CHPRCT

The Stainless Steel King 3 Tier Trolley Utility Cart is a versatile and durable addition to any home or professional setting.

3 shelves
 wheeled

Barcode: 0710579687341
 Nappi Code:N/A

Cart Stainless Steel - 3 Shelf



CH741DT

2 drawers and cupboard.
 Cardiac board.
 Infusion stand.
 oxygen bottle holder.
 Sliding top with drug compartments.
 Lower cupboard (one side only)

Barcode: 6009701458237
 Nappi Code:N/A

Medicine Trolley - Double Sided



CHPRCT

ABS Top Board
 Cold steel Spray Material
 Product Size:650x480x920mm

Barcode: 6009701459845
 Nappi Code:N/A

Procedure Cart - ABS Top Board

Trolleys



CHTT01

Stainless steel frame
4 wheels with brake
height adjustable
Gross Weight: 5kg

Barcode: 6009701459357
Nappi Code:N/A

Stainless Steel Tray Trolley



CHMRC

Medical Record Cart

Barcode: 0710579688553
Nappi Code: N/A

Medical Record Cart




CHM6D

Material: Stainless steel 304 1.0mm thickness
4 wheels with brakes
6 drawers
Product Size: 71.5x46x106cm
New Weight: 59kg
Gross Weight: 65kg

Barcode:6009701459142
Nappi Code:N/A

Medical Trolley Detachable - 6 Drawer



CH386

Detachable type
Material: stainless steel
2 shelves, 1 big drawer, 4 wheels with brakes
Product size: 66*44*86cm
Net weight: 12kg

Barcode:
Nappi Code:N/A

Stainless Steel Trolley Full Drawer - Anaesthetic



CHDC4

Lockable
PP material
4 shelves
Product size: 890*590*460mm
Net weight: 35kg

Barcode: 6009701459364
Nappi Code:N/A

Drug Cabinet - Lockable

Crutches and Walking Sticks



CH920

A walking stick is a device used primarily to aid walking, provide postural stability or support, or assist. Stability
Suitable for aged persons with poor balance
Height can be adjusted easily with a push button
Lightweight and sturdy
Good grip on any surface

- Total Height: 71 – 94cm
- Max Weight: 130kg

Barcode :6009701451825
Nappi Code:1103437001


Walking Stick Std



CHEB 12mm Ferrules

Barcode :
6009701457308

12mm Ferrules




CH933L(S)(M)(L)

Small: Total Height – 79 – 109cm
Medium: Total Height – 87 – 117cm
Large: Total Height – 95 – 125cm

Max Weight: 130kg

Barcode :
6009701452730
6009701452747
6009701452754
Nappi Code:
1103447001
1103450001
1103452001

Elbow Crutch-Small/ Med/ Large



CH935L(S)(M)(L)

- Semi Adjustable
- Lightweight wood construction
- Tough and durable
- Large non-slip rubber tips
- Maximum user weight of 150kg

Barcode :
6009701452532
6009701452501
6009701452556
Nappi Code:
1106515001
1106525001
1106538001

Wooden Crutches 1 Pair - Small / Medium / Large



CH925(S)(M)(L)

- Semi Adjustable
- Lightweight wood construction
- Tough and durable
- Large non-slip rubber tips
- Maximum user weight of 150kg

Barcode :
6009701453171
6009701453157
6009701453164
Nappi Code:
1103902001
1103902001
1103902001

UnderArm Crutch Aluminium Small/Med/Large 1 Pair



CH9274L

A foldable walking stick which makes it easy to pack away and uses up less space.

Max Weight 130kg

Total Height:83.5 – 93.5 cm

Barcode :6009701455540
Nappi Code:1103438001

Walking Stick Foldable with torch



CH926

- 3 legs at base of stick to give extra support.
- Comfortable grip handle
- Lightweight
- Max Weight: 110kg
- Total Height: 73 – 96cm

Barcode :6009701452464
Nappi Code:1034441001

Walking Stick Tripod



CH9207LW-3

- Walking Stick – Tripod
- Height Adjustable Walking Stick with 3 Legs for Elderly – Tripod
- 3 leg design for better stability
- Suitable for aged persons with poor balance
- Height can be adjusted easily with a push button
- Lightweight and sturdy

Barcode :6009701453553
Nappi Code:1103977001

Tripod Crutch

Crutches and Walking Sticks



CH9207LW-2

- Walking Stick – Tripod
- Height Adjustable Walking Stick with 3 Legs for Elderly – Tripod
 - 3 leg design for better stability
 - Suitable for aged persons with poor balance
 - Height can be adjusted easily with a push button
 - Lightweight and sturdy

Barcode:6009701455656
Nappi Code:1034441001

Walking Stick Tripod Compact

CH911L



- Folding and portable design.
Three feet support
Handle is made of polymer foam material, comfortable grip.
When used as a cane, the cane is 870mm / 34 tall.
The seat is approximately 500mm/20 high off the ground, perfect for comfortable sitting.
Non-slip feet, equipped with anti skid rubber floor mat, provide a good ability of the friction.

Barcode:6009701451924
Nappi Code:1106539001

Walking Stick with Seat



CH937-1

Elbow Crutch is a mobility aid that help to transfer the body's weight from an injured leg/foot to the upper body.
The user will slip their forearm into a cuff and then grasp a hand grip.

- Total Height: 95 – 118 cm
- Max Weight: 100kg

Barcode:6009701452372
Nappi Code:1103456001

Elbow Crutch Black Single



CH937-2

Elbow Crutch is a mobility aid that help to transfer the body's weight from an injured leg/foot to the upper body.
The user will slip their forearm into a cuff and then grasp a hand grip.

- Total Height: 95 – 118 cm
- Max Weight: 100kg

Barcode:6009701452389
Nappi Code:1103456001

Elbow Crutch Red Single



CH9314L

A walking stick with a sturdy base.
Rubber Grip comfortable handle with velcro strap for extra support

- Total Height:105 – 128cm
- Weight Capacity: 130kg

Barcode:6009701453676
Nappi Code:1103443001

Walking Stick with Velcro Arm Strap



CH936L

- New joint design provides superior tactile transmission.
- Cane fold and unfold with ease.
- Rolling tip, folds down 4 sections.
- Double elastic cord assembly for added security
- Hook style tip- with the security of the tip held in place with the elastic cord.

Barcode:6009701455625
Nappi Code:1120529001

White Cane - Blind Cane



CH9287L

A walking stick with a sturdy base.
Rubber Grip comfortable handle
4 Feet to support

- Total Height: 64 – 100cm
- Weight Capacity: 130kg

Barcode:6009701455663
Nappi Code:1103978001

Quadripod Walking stick with torch

Crutches and Walking Sticks



CH922 A walking stick with a sturdy base. Rubber Grip comfortable handle
4 Feet to support

- Total Height: 64 – 87cm
- Weight Capacity: 130kg

Barcode:6009701451092
Nappi Code:1103464001

Walking Stick Quadripod Grey



CH924 A walking stick with a sturdy base. Rubber Grip comfortable handle
4 Feet to support

- Total Height: 64 – 100cm
- Weight Capacity: 130kg

Barcode:6009701451917
Nappi Code:1118679001

Walking Stick Quadripod -Extra Length



CH924L A walking stick with a sturdy base. Rubber Grip comfortable handle
4 Feet to support

- Total Height: 64 – 87cm
- Weight Capacity: 130kg

Barcode:6009701451917
Nappi Code:1118679001

Walking Stick Quadripod -Extra Length Black



CH924F Adjustable Metal Folding Cane
Comfortable contour grip with rubber tip
Height adjusts from 33inch to 37inches
in 1 inch increments
4 folding sections snap out
automatically

Anodized aluminum shaft

Barcode:0710579687358
Nappi Code:N/A

Walking Stick Std Foldable



CH937-3

Elbow Crutch is a mobility aid that help to transfer the body's weight from an injured leg/foot to the upper body. The user will slip their forearm into a cuff and then grasp a hand grip.

- Total Height: 95 – 118 cm
- Max Weight: 100kg

Barcode:
Nappi Code:

Elbow Crutch Orange Single



CLINI *health*
leaders in healthcare



ORTHOPAEDICS



Knee and Leg Support



Neoprene Hinge Knee Support allows compression around the knee, with the support of side hinges it allows normal flexion of the knee. Extra strap allows secured fitting and hinges mimics the motion of a natural knee joint.

Fibre Contents: Polyester, Polypropylene, Rubber & Rubber

Size: Small
Product Code: NP050RCHS
Barcode:
Nappi Code:

Size: Medium
Product Code: NP050RCHM
Barcode:
Nappi Code:

Size: Large
Product Code: NP050RCHL
Barcode:
Nappi Code:

Size: XL
Product Code: NP050RCHXL
Barcode:
Nappi Code:

Size: XXL
Product Code: NP050RCHXXL
Barcode:
Nappi Code:

Knee Support Hinge Neoprene - S-XXL

High-performance design for active support. Dual strapping provides customized compression. Soft material for all-day comfort. Sleek design fits inside the shoe.



Barcode: 8906164422905
Nappi Code: 1175245001

Ankle Wrap - Children Universal

High-performance design for active support. Dual strapping provides customized compression. Soft material for all-day comfort. Sleek design fits inside the shoe.



Barcode: 8906164422905
Nappi Code: 1175246001

Ankle Wrap Neoprene



Knee Immobilize is specially designed to provide rigid support and protection to the injured knee, mostly post-operative cases. Extra Patellar strapping enhances the immobilization of the knee.

FIBRE CONTENTS: Nylon, Cotton, Polyester & High Density Foam

Size: Small
Product Code: KS010RCHS
Barcode: 8906164421397
Nappi Code: 1171652001

Size: Medium
Product Code: KS010RCHM
Barcode: 8906164421403
Nappi Code: 1171655001

Size: Large
Product Code: KS010RCHL
Barcode: 8906164421410
Nappi Code: 1171658001

Size: XXLarge
Product Code: KS010RCHXXL
Barcode: 8906164421434
Nappi Code: 1171663001

Size: X-Large
Product Code: KS010RCHXL
Barcode: 8906164421427
Nappi Code: 1171661001

Knee Immobiliser 22inch



Knee Immobilize is specially designed to provide rigid support and protection to the injured knee, mostly post-operative cases. Extra Patellar strapping enhances the immobilization of the knee.

FIBRE CONTENTS: Nylon, Cotton, Polyester & High Density Foam

Barcode: 8906164422516
Nappi Code: 1171732001



KS030RCHC

Knee Immobiliser - Child

Knee and Leg Support



Knee Immobilize is specially designed to provide rigid support and protection to the injured knee, mostly post-operative cases. Extra Patellar strapping enhances the immobilization of the knee.

FIBRE CONTENTS: Nylon, Cotton, Polyester & High Density Foam

Size: Small
Product Code: KS02ORCHS
Barcode: 8906164421458
Nappi Code: 1171636001

Size: Medium
Product Code: KS02ORCHM
Barcode: 8906164421465
Nappi Code: 1171639001

Size: Large
Product Code: KS02ORCHL
Barcode: 8906164421472
Nappi Code: 1171642001

Size: XXL Large
Product Code: KS02ORCHXXL
Barcode: 8906164421496
Nappi Code: 1171648001

Size: X-Large
Product Code: KS02ORCHXL
Barcode: 8906164421489
Nappi Code: 1171646001

Knee Immobiliser 19inch

ROM Knee Brace controls and restricts the degree flexion of the knee joint through multi-point fixation system to allow a range of motion of immobilization to the knee. ROM dial with 0 - 120 degree flexion and 0 - 60 degree hyper extension, provides locking at any angle or motion between any two angles with intervals of 15 degrees.

FIBRE CONTENTS: Polyester, Foam

Barcode: 8906164421885
Nappi Code: 1171724001

KS11ORCHU



ROM Knee Brace - Universal

Neoprene Knee Wrap is an open design product which offers easy application on swollen or asymmetric knee and allows customized compression. Neoprene fabric provides firm compression, warmth & support to the knee to allay pain and inflammation.

Fibre Contents: Polyester, Rubber, Polypropylene & Nylon

Barcode: 8906164422479
Nappi Code: 1172565001



NP04ORCHU

Knee Wrap Neoprene - Universal



NP04CORCHU

Neoprene Knee Wrap is an open design product which offers easy application on swollen or asymmetric knee and allows customized compression. Neoprene fabric provides firm compression, warmth & support to the knee to allay pain and inflammation.

Fibre Contents: Polyester, Rubber, Polypropylene & Nylon

Barcode: 8906164422479
Nappi Code: 1172566001



Knee Wrap - Children Universal



CLINIhealth
leaders in healthcare

Knee and Leg Support



FIBRE CONTENTS: Nylon, Cotton, Polyester & High Density Foam

Size: Small
Product Code: KS13ORCHS
Barcode:
Nappi Code:

Size: Medium
Product Code: KS13ORCHM
Barcode:
Nappi Code:

Size: Large
Product Code: KS13ORCHL
Barcode:
Nappi Code:

Size: X-Large
Product Code: KS13ORCHXL
Barcode:
Nappi Code:

Hinge Knee Brace



Size: Small
Product Code: KS16ORCHS
Barcode:
Nappi Code:

Size: Medium
Product Code: KS16ORCHM
Barcode:
Nappi Code:

Size: Large
Product Code: KS16ORCHL
Barcode:
Nappi Code:

Size: XL
Product Code: KS16ORCHXL
Barcode:
Nappi Code:

Size: XXL
Product Code: KS16ORCHXXL
Barcode:
Nappi Code:

Tubular Knee Support Beige



Barcode:
Nappi Code:

KS12ORCHU

Patellar Support - Universal

Knee and Leg Support

Tubular Knee Support helps provide extra support & warmth compression to the knee joint. Tubular design helps you to wear it whole day without feeling any discomfort

Size: Small
Product Code: KS08ORCHS
Barcode:8906164421663
Nappi Code: 11717380011

Size: Medium
Product Code: KS08ORCHM
Barcode:8906164421670
Nappi Code: 1171740001

Size: Large
Product Code: KS08ORCHL
Barcode:8906164421687
Nappi Code: 1171742001

Size: X Large
Product Code: KS08ORCHXL
Barcode:8906164421694
Nappi Code: 1171744001

Size: XXL Large
Product Code: KS08ORCHXXL
Barcode:8906164421700
Nappi Code: 1171746001



Tubular Knee Support Open Patella

Back and Shoulder Support

Contoured Lumbo Sacral Support is designed to provide support & comfort to the patient in low back pain. With contoured aluminum splints, it provides support to weak back muscles & corrects the postural deformity.

FIBRE CONTENTS: Polyester & Rubber

Size: Small
Product Code: BS04ORCHS
Barcode:8906164420406
Nappi Code: 1171007001

Size: Medium
Product Code: BS04ORCHM
Barcode:8906164420413
Nappi Code: 1171009001

Size: Large
Product Code: BS04ORCHL
Barcode:8906164420420
Nappi Code: 1171011001

Size: XL
Product Code: BS04ORCHXL
Barcode:8906164420437
Nappi Code: 1171014001

Size: XXL
Product Code: BS04ORCHXXL
Barcode:8906164420444
Nappi Code: 1171015001



Contoured LS Support Belt



Posture Corrector is an effective device to maintain correct posture during daily activities and long sitting by applying a correctional force. Posture Corrector belt has customizable design and elasticity of materials to ensure correctional force and higher comfort and for effective use of this product, ensure a strong, snug and comfortable fit.

Size: Small
Product Code: BS17ORCHS
Barcode:8906164420680
Nappi Code: 1171217001

Size: Medium
Product Code: BS17ORCHM
Barcode:8906164420697
Nappi Code: 1171221001

Size: Large
Product Code: BS17ORCHL
Barcode:8906164420703
Nappi Code: 1171224001

Size: XL
Product Code: BS17ORCHXL
Barcode:8906164420710
Nappi Code: 1171226001

Size: XXL
Product Code: BS17ORCHXXL
Barcode:8906164420727
Nappi Code: 1171227001

Posture Corrector

Back and Shoulder Support



Rigid Splints are light in weight, durable and provide effective immobilization. They can be customized to provide accurate shape and excellent comfort. They maintain the spine in neutral position and allow accurate posture correction.

Horizontal Chest Splints arrest the rotational movement of the torso for effective immobilization.

Size: Universal
Product Code: FA07ORCHU
Barcode: 8906164420994
Nappi Code: 1171713001

Size: Special
Product Code: FA07ORCHS
Barcode: 8906164421007
Nappi Code: 1171716001



Taylor Brace

Soft PU foam padding ensures no underarm rash, improves comfort and compliance. Figure of eight design with lateral vector ensures linear union of the clavicle bone, allows the required movement of the arms with perfect immobilization of the clavicle bone, easy application and tightening of the product.

Size: Medium
Product Code: FA04ORCHM
Barcode: 8906164420888
Nappi Code: 1171821001

Size: Large
Product Code: FA04ORCHL
Barcode: 8906164420895
Nappi Code: 1171822001



Clavicle Brace



Anatomic design maintains the natural low back curve and relieves the anterior pressure on the vertebrae. Side supports hold the body and enhance the comfort.

Barcode: 8906164420529
Nappi Code: 1171824001



Back Rest - Universal

Elastic Shoulder Immobilizer is a lightweight, sleek and smart design to provide immobilization in shoulder dislocation and post-operative rehabilitation of the shoulder. It reduces abduction and arm rotation by positioning the arm close to the body.

FA05ORCHL



Barcode: 8906164420932
Nappi Code: 1171718001

Elastic Shoulder Immobiliser - Large


Two-piece design is more comfortable, offers better immobilization, easy application and removal. Can be tightened as per need for enhanced immobilization and grip.

Size: Universal

Size: Child

Product Code: FA08ORCHU
Barcode: 8906164421014
Nappi Code: 1171719001

Product Code: FA08ORCHC
Barcode: 8906164422592
Nappi Code: 1171720001



Universal Shoulder Immobiliser

Back and Shoulder Support

Anatomical shape ensures good grip around the shoulder, perfect fitting and optimal immobilization. Stretchable shoulder band ensures strong grip on the shoulder.

Size: Medium
Product Code: NP06ORCHM
Barcode: 8906164422844
Nappi Code:

Size:XL
Product Code: NP06ORCHXL
Barcode: 8906164422868
Nappi Code:

Size:Large
Product Code: NP06ORCHL
Barcode: 8906164422851
Nappi Code:

Size:XXL
Product Code: NP06ORCHXXL
Barcode: 8906164422875
Nappi Code:



Shoulder Support

Arm, Wrist, Hand and Foot Support



Recuperate in comfort with our Arm Sling. It provides you a comfortable protection & relief to fractured arm while you recuperating & healing. Simple closure with straps helps you to wear it easily.

FIBRE CONTENTS: Cotton, Polyester & Medium Density Foam

Size: Small
Product Code: FA01ORCHS
Barcode: 8906164420789
Nappi Code: 1170995001

Size: XL
Product Code: FA01ORCHXL
Barcode: 8906164420819
Nappi Code: 1171001001

Size: Medium
Product Code: FA01ORCHM
Barcode: 8906164420796
Nappi Code: 1170997001

Size: Child
Product Code: FA01ORCHCH
Barcode: 8906164422523
Nappi Code: 1171002001

Size: Large
Product Code: FA01ORCHL
Barcode: 8906164420802
Nappi Code: 1170999001



Arm Sling

Rocker sole deloads the ankle joint, enhances the ease of walking on foot cast, restores the normal gait. It improves balance of patient on foot cast and helps in early ambulation. It is extremely light in weight, colorful and aesthetically pleasing.

Size: Small
Product Code: FA03ORCHS
Barcode: 8906164420833
Nappi Code: 1171352001

Size: XL
Product Code: FA03ORCHXL
Barcode: 8906164420864
Nappi Code: 1171358001

Size: Medium
Product Code: FA03ORCHM
Barcode: 8906164420840
Nappi Code: 1171355001

Size: Large
Product Code: FA03ORCHL
Barcode: 8906164420857
Nappi Code: 1171357001



Cast Shoe

Arm, Wrist, Hand and Foot Support

Removable, anatomically shaped, double splinting ensures rigid immobilization. It can be customized to provide the required degree of wrist dorsiflexion. Lace pull mechanism allows tightening and loosening of the brace with a single pull stroke. Easy to wear

Size: Universal

Size: Child

Product Code: WS100RHCU
Barcode: 8906164421274
Nappi Code: 1171721001

Product Code: WS110RCH
Barcode: 8906164422660
Nappi Code: 1171722001



Wrist and Forearm Brace

Designed with a dorsal shell made of polypropylene with a leaf spring action to arrest the plantar flexion movement and lift the foot.

Size: Left Medium
Product Code: ASL04ORCHM
Barcode: 8906164422042
Nappi Code:

Size: Left Large
Product Code: ASL04ORCHL
Barcode: 8906164422059
Nappi Code:

Size: Right Medium
Product Code: ASR04ORCHM
Barcode: 8906164422080
Nappi Code:

Size: Right Large
Product Code: ASR04ORCHL
Barcode: 8906164422097
Nappi Code:



Foot Drop Splint

NP01ORCHU



Removable, anatomically shaped, double splinting ensures rigid immobilization. It can be customized to provide the required degree of wrist dorsiflexion. Lace pull mechanism allows tightening and loosening of the brace with a single pull stroke. Easy to wear and remove. Reverse buckle mechanism ensures controlled tightening and better immobilization. Ergonomic and sleek design.

Barcode: 8906164422448
Nappi Code:



Wrist Binder Neoprene - Universal

Arm, Wrist, Hand and Foot Support



Diabetic Socks

GP15ORCHM - Medium
GP15ORCHL - Large

Barcode: 8906164423179
Nappi Code: 1172573001

Diabetic Socks - Universal

Customizable Splints.
Anatomical thumb abduction.
Fits both left and right hand.

NP01ORCHU



Barcode: 8906164420932
Nappi Code: 1171718001

Wrist Binder Thumb Support - Universal

Neck & Coccyx Supports



Cervical Collar provides a comfortable & warm support to your neck pain & stiffness. Its skin friendly fabric with sweat absorbing quality gives you a pleasant feel while healing your pain at a faster pace.

FIBRE CONTENTS: Cotton, Polyester & Medium Density Foam


<p>Size: Small Product Code: CS02ORCHS Barcode: 8906164420055 Nappi Code:</p>	<p>Size: Medium Product Code: CS02ORCHM Barcode: 8906164420062 Nappi Code:</p>	<p>Size: Large Product Code: CS02ORCHL Barcode: 8906164420079 Nappi Code:</p>	<p>Size: X-Large Product Code: CS02ORCHXL Barcode: 8906164420086 Nappi Code:</p>
--	---	--	---

Cervical Collar Soft

Soft polyurethane body provides unmatched comfort. It is used as prophylactic collar, reduces pain and provides comfortable support during traveling.

Covered with hypoallergenic cotton stockinette improves comfort.

<p>Size: Small Product Code: CS05ORCHS Barcode: 8906164420147 Nappi Code: 1171038001</p>	<p>Size: Medium Product Code: CS05ORCHM Barcode: 8906164420154 Nappi Code: 1171042001</p>	<p>Size: Large Product Code: CS05ORCHL Barcode: 8906164420161 Nappi Code: 1171044001</p>
<p>Size: X-Large Product Code: CS05ORCHXL Barcode: 8906164420178 Nappi Code: 1171045001</p>	<p>Size: Child Product Code: CS05ORCHC Barcode: 8906164422493 Nappi Code:</p>	



Boneless Collar Soft

Neck & Coccyx Supports



Coccyx Cushion is anatomically to provide maximum comfort while reducing pressure on the tailbone. It is made up with PU foam that helps in attaining the required support while sitting for long duration. It is ideal for pregnant women, people with spondylitis, sacral joint pain, sciatica or stenosis.

FIBRE CONTENTS: Polyester, PU Foam, Lycra

Barcode: 8906164420536
Nappi Code: 1171825001



BS09ORCHU

Coccyx Cushion - Universal

Cervical Pillow is designed to keep the neck in the slightly hyper extension position which provides optimal muscle relaxation and proper support to head and neck during sleep.



CS06ORCHU

Barcode: 8906164420185
Nappi Code:



Contoured Cervical Pillow - Universal

Hernia Belts

Offers a gentle and constant relief from reducible hernia by providing adjustable and comfortable pressure to the hernia.

On wearing the belt the occurrence of herniation gets eliminated.

Hernia Belt is available in skin color, and can be worn inconspicuously under clothing.
Molded foam Pads
Adjustable leg and abdominal strap
Soft and skin friendly materials

Size: Medium
Product Code: BS06ORCHM
Barcode: 8906164420482
Nappi Code: 1171726001

Size: Large
Product Code: BS06ORCHL
Barcode: 8906164420154
Nappi Code: 1171728001

Size: XL
Product Code: BS06ORCHXL
Barcode: 8906164420505
Nappi Code: 1171730001



Hernia Belt

Compression Stockings

Easy Conform Material
Helps Relieve Varicose Veins
Graduation Compression
Improves Circulation

Size: Small

Product Code: COMPSTC2BKSL
Barcode: 6009701458961

Size: Medium

Product Code: COMPSTC2BKM
Barcode: 6009701458978

Size: Large

Product Code: COMPSTC2BKL
Barcode: 6009701458985

Size: XL


Product Code: COMPSTC2BKXL
Barcode: 6009701458992

Size: XXL

Product Code: COMPSTC2BKXXL
Barcode: 6009701459043

Size: XXXL

Product Code: COMPSTC2BKXXXL
Barcode: 6009701455588



COMPSTC2

Compression Stockings Below the Knee Class 2

Easy Conform Material
Helps Relieve Varicose Veins
Graduation Compression
Improves Circulation

Size: Small

Product Code: COMPSTC3BKSL
Barcode: 6009701459005

Size: Medium

Product Code: COMPSTC3BKM
Barcode: 6009701459012

Size: Large

Product Code: COMPSTC3BKL
Barcode: 6009701459029

Size: XL

Product Code: COMPSTC3BKXL
Barcode: 6009701459036

Size: XXL

Product Code: COMPSTC3BKXXL
Barcode: 6009701459050

Size: XXXL

Product Code: COMPSTC3BKXXXL
Barcode: 6009701455595



COMPSTC3

Compression Stockings Below the Knee Class 3

Foot and Ankle



Application of this device is recommended only when the fracture is stable and there are acceptable limits of angular and rotational deformity.

The frequency and duration of the use should be determined by your prescribing health care professional.

Size: Small

Product Code: AIRBTLGCH
Barcode: 6009701458015
Nappi Code: 1156306001

Size: Medium

Product Code: AIRBTLMGCH
Barcode: 6009701458008
Nappi Code: 1156307001

Size: Large

Product Code: AIRBTLGCH
Barcode: 6009701457988
Nappi Code: 1156308001

Size: X-Large

Product Code: AIRBTLXGCH
Barcode: 6009701457971
Nappi Code: 1156609001

Air Boot Long Adjustable Fixation Brace DH Walker VIII Grey

Application of this device is recommended only when the fracture is stable and there are acceptable limits of angular and rotational deformity.

The frequency and duration of the use should be determined by your prescribing health care professional.

Size: Small

Product Code: AIRBTLBCH
Barcode: 6009701458107
Nappi Code: N/A

Size: Medium

Product Code: AIRBTLMBCH
Barcode: 6009701458114
Nappi Code: N/A

Size: Large

Product Code: AIRBTLBCH
Barcode: 6009701458121
Nappi Code: N/A

Size: X-Large

Product Code: AIRBTLXLBCH
Barcode: 6009701458138
Nappi Code: N/A



Air Boot Long Adjustable Fixation Brace DH Walker VIII Black



Foot and Ankle

Application of this device is recommended only when the fracture is stable and there are acceptable limits of angular and rotational deformity.

The frequency and duration of the use should be determined by your prescribing health care professional.

Size: Small Product Code: ADJANFSCH Barcode: 6009701457551 Nappi Code: 1156336001	Size: Medium Product Code: ADJANFMCH Barcode: 6009701457568 Nappi Code: 1156337001	Size: Large Product Code: ADJANFLCH Barcode: 6009701457575 Nappi Code: 1156338001
Size: XL Product Code: ADJANFXLCH Barcode: 6009701457582 Nappi Code: 1156339001		



Walker Adjustable Ankle and Foot Fixation DH Walker Long Type with Airbag

Providing stabilizing support following an injury or surgery, the knee brace features a strong carbon fiber splint and adjustable bilateral range of motion settings to prevent re-injury and help facilitate rehabilitation. The reversible knee brace allows the flexion and extension to be set to limit motion or support a full range of motion for active individuals.

Size: Left Medium Product Code: WLADJLMCH Barcode: 6009701457629 Nappi Code: 1156515001	Size: Left Large Product Code: WLADJLLCH Barcode: 6009701457636 Nappi Code: 1156521001
Size: Right Medium Product Code: WLADJRMCH Barcode: 6009701457643 Barcode: 1156523001	Size: Right Large Product Code: WLADJRLCH Barcode: 6009701457650 Nappi Code: 1156526001



Walker Adjustable Ankle and Foot Brace High Model IV



Designed for rehabilitation after injury, fracture, sprains or surgery of foot, ankle or lower leg to prevent plantar flexion, when required. The rigid casing provides support to lift the ankle, constraining the plantar flexion, enabling prevention of conditions like foot drop or contractures on prolonged bed rest.

Size: Small Product Code: WADJADSCH Barcode: 6009701458046 Nappi Code: N/A	Size: Medium Product Code: WADJADMCH Barcode: 6009701458053 Nappi Code: N/A	Size: Large Product Code: WADJADLCH Barcode: 6009701458060 Nappi Code: N/A
--	---	--

Walker Adjustable Ankle and Foot Brace High Model Night Splint III

Designed for rehabilitation after injury, fracture, sprains or surgery of foot, ankle or lower leg to prevent plantar flexion, when required. The rigid casing provides support to lift the ankle, constraining the plantar flexion, enabling prevention of conditions like foot drop or contractures on prolonged bed rest. De-rotation Bar can be used for prevention of rotation post femur fracture.

Size: Left Small Product Code: WAADJCHLSCH Barcode: 6009701457667 Nappi Code: 1156532001	Size: Left Medium Product Code: WAADJCHLMCH Barcode: 6009701457674 Nappi Code: 1156533001	Size: Left Large Product Code: WAADJCHLLCH Barcode: 6009701457001 Nappi Code: N/A
Size: Right Small Product Code: WAADJCHRSCCH Barcode: 6009701458077 Barcode: N/A	Size: Right Medium Product Code: WAADJCHRMCH Barcode: 6009701458084 Nappi Code: N/A	Size: Right Large Product Code: WAADJCHRLCH Barcode: 6009701458091 Nappi Code: N/A



Walker Adjustable Ankle and Foot Brace Low Model- Child

Foot and Ankle



Intended to prevent or treat acute ankle injuries, chronic instability, post-operative use and ankle sprains grade I, II, and III
Semi-rigid, anatomically designed shells are lined with patented Duplex air cells, which help stabilize and protect the ankle

ANBRCACH



Barcode: 6009701457681
Nappi Code: 1156517001

Ankle Brace with clips and Airbag - Universal

Leg



Controls and restricts degree of flexion of the knee joint through a multi-point fixation system to allow a range of motion or immobilization to the knee.

Universal design allows it to be used on most patients and for multiple orthopaedic problems.

Colour: Black
Product Code : ROMKN120BCH
Barcode: 6009701450491
Nappi Code: 1156522001

Colour: Grey
Product Code : ROMKN120GCH
Barcode: 6009701457711
Nappi Code: 1156522001



ROM Knee Brace Fixed Support 10 - 120 degrees II Grey / Black



Controls and restricts degree of flexion of the knee joint through a multi-point fixation system to allow a range of motion or immobilization to the knee.

Universal design allows it to be used on most patients and for multiple orthopaedic problems.

Barcode: 6009701457728
Nappi Code: 1156525001



ROMKN90CH

ROM Knee Brace Fixed Support 10 - 90 degrees V - Universal



Open design, easy to adjust the size and elastic,
Use strong, durable material imported from Japan, feel soft and comfortable,
Easy to adjust the knee joint of aluminum alloy stent, can very good for control of knee flexion angles.

Barcode: 6009701452907
Nappi Code: N/A

ADJKFXLCH

Adjustable Knee Brace Fixed Support XII Large



Back

High Quality Durable Women Maternity Support Belt Pregnancy Belly Back Brace Belt
Pregnancy belt supports the abdominal area and helps ease the discomfort of pregnancy by supporting your baby bump.

Size: Small
Product Code: MATSPBSCH
Barcode: 6009701457735
Nappi Code: 1156520001

Size: Medium
Product Code: MATSPBMCH
Barcode: 6009701457742
Nappi Code: 1156524001

Size: Large
Product Code: MATSPBLCH
Barcode: 6009701457759
Nappi Code: 1156527001

Size: XL
Product Code: MATSPBXLCH
Barcode: 6009701457766
Nappi Code: 1156529001



Maternity Support Belt

Hand

Thumb Stabilization and Pain Relief- This thumb brace can help with arthritis and other chronic conditions, but is also well suited for sprains, fractures, or even soft tissue injuries to the thumb

Fingers Are Free in Motion- stabilizes the thumb while giving free range of movement to all the other digits, like carrying a child, typing, playing the instruments, doing sports, cargo handling etc.



Size: Left Small
Product Code: THBRSSLCH
Barcode: 6009701457810
Nappi Code: 1156443001

Size: Left Medium
Product Code: THBRSMMLCH
Barcode: 6009701457827
Nappi Code: 1156445001

Size: Left Large
Product Code: THBRSLLLCH
Barcode: 6009701457834
Nappi Code: 1156448001

Size: Right Small
Product Code: THBRSSRCH
Barcode: 6009701457841
Barcode: 1156450001

Size: Right Medium
Product Code: THBRSMRCH
Barcode: 6009701457858
Nappi Code: 1156452001

Size: Right Large
Product Code: THBRSLLRH
Barcode: 6009701457865
Nappi Code: 1156453001

Thumb Brace with Strap

Breathable mesh fabric for ultimate comfort
Adjustable straps provide customized fit and optimal pressure
Built-in aluminum sheet provides strong support
Suitable for both left and right hands

Size: Left Small
Product Code: CAFXSTSLCH
Barcode: 6009701457599
Nappi Code: 1156535001

Size: Left Medium
Product Code: CAFXSTMLCH
Barcode: 6009701457605
Nappi Code: 1156537001

Size: Left Large
Product Code: CAFXSTLLCH
Barcode: 6009701457612
Nappi Code: 1156539001

Size: Right Small
Product Code: CAFXSTSRCH
Barcode: 6009701458145
Barcode: N/A

Size: Right Medium
Product Code: CAFXSTMTRCH
Barcode: 6009701458152
Nappi Code: N/A

Size: Right Large
Product Code: CAFXSTLRCH
Barcode: 6009701458169
Nappi Code: N/A



Carpal Fixation Strap

Hand

Finger splints are used to keep these vulnerable extremities in place while you heal from a finger sprain or break, or to help increase a limited range of motion from a chronic condition. Aids bruised, sprained or cut finger by immobilizing movement to protect the injured finger. Features and Benefits. Designed to help relieve finger pain and reduce the swelling caused by sprains, jams, strains and arthritis.

Size: Small
Product Code: FINSPSCH
Barcode: 6009701457872
Nappi Code: 1156488001

Size: Medium
Product Code: FINSPMCH
Barcode: 6009701457889
Nappi Code: 1156489001

Size: Large
Product Code: FINSPLCH
Barcode: 6009701457896
Nappi Code: 1156490001



Finger Splint

Arm



Controls and restricts degree of flexion and extension of the elbow joint through sophisticated dial system to allow a range of motion or immobilization of the elbow. Universal design allows it to be used on most of the patients and for multiple orthopaedic problems.

Adjustable arm length
Universal size

Size: Left
Product Code: ADELBRLCH
Barcode: 6009701457698
Nappi Code: 1156518001

Size: Right
Product Code: ADELBRRCH
Barcode: 6009701457704
Nappi Code: 1156519001

Adjustable Elbow Brace

Arm Sling hold, support and the sprained, broken or surgically operated arm in the flexion position while it recuperates.

Designed to cover the cast and give extra protection to the injured arm.

Adjustable

Enhanced protection

Size: Small
Product Code: ARSLINSCH
Barcode: 6009701457780
Nappi Code: 1156499001

Size: Medium
Product Code: ARSLINMCH
Barcode: 6009701457797
Nappi Code: 1156500001

Size: Large
Product Code: ARSLINLCH
Barcode: 6009701457803
Nappi Code: 1156501001



Arm Sling



CLINIhealth
leaders in healthcare

Shoulder

Application: shoulder pain, joint weakness, minor contusion/sprain, joint dislocation/instability.

Features:

1. Adjustable design, suitable for most people
2. Using seamless sewing technology
3. Breathable
4. Comfortable to wear

Size: Left

Product Code: SHBRAILCH

Barcode: 6009701457902

Nappi Code: 1156491001

Size: Right

Product Code: SHBRAIRCH

Barcode: 6009701457919

Nappi Code: 1156492001



Shoulder Brace with Airbag







CLINI *health*®

 **IS AN ISO 13485
CERTIFIED COMPANY**



CONTACT DETAILS

Calibra Corner No 26 Michelson, Cnr Michelson & Kirschner Rd, Westwood, Boksburg, Unit D3

011 421 0702 • 011421 6989 • 084 245 7866 • 074 772 4132 • 084 366 8958

admin@clinihealth.co.za

www.clinihealth.co.za